

## One-step identification of five prominent chicken

### *Salmonella* serovars and biotypes

Chunhong Zhu<sup>1,†\*</sup>, Min Yue<sup>1\*</sup>, Shelley Rankin<sup>1</sup>, François-Xavier Weill<sup>2</sup>, Joachim Frey<sup>3</sup>, and Dieter M. Schifferli<sup>1#</sup>

<sup>1</sup>Department of Pathobiology, University of Pennsylvania School of Veterinary Medicine, Philadelphia, PA, USA

<sup>2</sup>Institut Pasteur, Unité des Bactéries Pathogènes Entériques, Paris, France

<sup>3</sup>Institute of Veterinary Bacteriology, Department of Infectious Diseases and Pathobiology, Vetsuisse Faculty, University of Bern, Bern, Switzerland.

<sup>#</sup>Correspondence should be addressed to D.M.S. ([dmschiff@vet.upenn.edu](mailto:dmschiff@vet.upenn.edu))

<sup>†</sup>Current address: Jiangsu Institute of Poultry Science, Chinese Academy of Agricultural Sciences, Yangzhou, 225125 Jiangsu, China

\*Chunhong Zhu and Min Yue contributed equally to this work.

#### Abstract:

Based on bacterial genomic data, we developed a one-step multiplex PCR assay to identify *Salmonella* and simultaneously differentiate the two invasive avian-adapted *S. enterica* serovar Gallinarum biotypes Gallinarum and Pullorum, as well as the most frequent, specific and asymptomatic colonizers of chickens, serovars Enteritidis, Heidelberg and Kentucky.

#### Keywords:

Chicken, *Salmonella*, *S. Enteritidis*, *S. Heidelberg*, *S. Kentucky*, *S. Gallinarum*, *S. Pullorum*, multiplex PCR

39 Strains of most *Salmonella* serovars are zoonotic. Approximately 90% of human  
40 salmonellosis results from ingestion of contaminated food products of animal or plant  
41 origin (1). With over 19,000 reported cases in the US for 2013, *Salmonella* remains  
42 the most frequently isolated bacterial food pathogen, as determined by the  
43 surveillance network FoodNet which pools the data of 10 US monitoring sites (2). In  
44 parallel to the rise of poultry consumption over the years in the US, the commercial  
45 poultry industry has grown impressively, reaching over 9 billion raised and processed  
46 broilers per year and a yearly production of over 77 billion table eggs, as indicated for  
47 2009 (3). *Salmonella* is a frequent asymptomatic intestinal colonizer of poultry. Stress  
48 or underlying diseases in young birds create optimal conditions for productive  
49 horizontal transmission of *Salmonella* sp. Data from the USDA-FSIS suggests that  
50 every fourth raw chicken part is likely contaminated with *Salmonella* (2). Moreover,  
51 major *Salmonella* serovars can spread to reproductive organs, leading to vertical  
52 transfer of the bacteria and egg-related salmonellosis (4, 5). Accordingly, poultry and  
53 egg consumption represent a significant source of *Salmonella* infections in the US.

54  
55 Four *Salmonella* serovars are of particular concern to the poultry industry, namely  
56 Enteritidis, Heidelberg, Kentucky and Gallinarum (6). *S. Gallinarum* is an invasive  
57 agent of chicken salmonellosis resulting in high mortality and morbidity, with biotype  
58 Pullorum (*S. Pullorum*) which causes "white diarrhea" in young chicken (pullorum  
59 disease), and biotype Gallinarum which is responsible for fowl typhoid (7). Although  
60 this serovar remains endemic in many countries, it has essentially been eradicated  
61 through culling programs in the domestic fowl industry of the USA and several other  
62 developed countries. *S. Gallinarum* can colonize and/or cause disease in various  
63 domestic and wild birds, which might explain its occasional detection in backyard  
64 birds of developed countries (8). In recent years, *S. Enteritidis* became a most  
65 frequently isolated serovar in poultry and from foodborne outbreaks linked to poultry  
66 products in developed countries (9). This serovar was suggested to have filled the  
67 ecological niche vacated by the eradicated *S. Gallinarum* biotypes Pullorum and  
68 Gallinarum (10). Lately, *S. Heidelberg* has become another major serovar responsible  
69 for foodborne infections from poultry products (11, 12), as well as one of the most  
70 common serovars obtained from non-clinical chicken isolates (9, 13, 14). *S. Kentucky*  
71 is the most common serotype isolated from chickens and the second most common  
72 one found among retail chicken product in the USA. It has been rarely reported in  
73 human cases in North America (15, 16), although this could change with worldwide  
74 spreading of the ciprofloxacin-resistant ST198 (17).

75  
76 Here we describe a simple one-step multiplex polymerase chain reaction (PCR)  
77 method to identify major chicken *S. enterica* subsp. *enterica* serovars. The approach  
78 was based on designing primers that specifically amplify unique sets of *Salmonella*  
79 spp. and serovar-associated DNA sequences in one PCR tube (Table 1), taking  
80 advantage of 3,161 available *Salmonella* genomes, including strains from serovar  
81 Enteritidis (369 genomes), Heidelberg (154), Kentucky (63), Gallinarum (8 biotype  
82 Pullorum and 4 biotype Gallinarum), and 2,563 genomes from 104 other serovars.

83 The desired specificities were checked by using BLAST (NCBI, non-redundant  
84 nucleotide collection). Strains of the *Salmonella* genus and Gallinarum biotypes were  
85 identified by primers for differently conserved DNA segments in the their *bcf* and *ste*  
86 fimbrial usher genes, respectively (18). Specific primers for serovar Gallinarum  
87 biotype Gallinarum were made by taking advantage of a deletion of 4 nucleotides in  
88 *steB* of biotype Pullorum. Other specific DNA signatures served as primer targets to  
89 separate serovars Enteritidis, Heidelberg and Kentucky. Briefly, for the multiplex PCR,  
90 pure template DNA (1-5 ng per reaction; MagNA Pure LC DNA Isolation Kit III,  
91 Roche Life Sci., Indianapolis, IN) or crude DNA (approximately 75 ng per reaction,  
92 from bacterial suspensions boiled for 5 min,  $10^7$  CFU/ $\mu$ l dH<sub>2</sub>O, using 1  $\mu$ l  
93 supernatants after centrifugation) was amplified with Taq DNA polymerase and a final  
94 concentration of 1.5 mM Mg<sup>2+</sup> (Choice Taq Blue™, Denville Sci. Inc., South  
95 Plainfield, NJ) using standard protocols. The PCR (25 cycles with an annealing  
96 temperature of 56°C) was performed with a Hybaid Thermal cycler (Thermo Fisher  
97 Sci., Waltham, MA). The specificity and compatibility of the primer sets in a  
98 multiplex PCR was assessed using genomic DNA from 128 *Salmonella* strains that  
99 included a total of 34 different serovars as well as three *Escherichia coli* and two  
100 *Yersinia* spp. as negative control strains (Suppl. Table 1).

101  
102 Agarose gel electrophoresis profiles for each different amplicon sets are visualized  
103 with representative strains in Fig. 1 and the results for all the strains are listed in Table  
104 2. All the *Salmonella* strains were recognized as such, as were strains of the  
105 Gallinarum biotypes and the Enteritidis, Heidelberg and Kentucky serovars. Thus, the  
106 obtained experimental results were in agreement with the genomic information used  
107 for the primers' design and validated the proposed identification of *S. enterica* and the  
108 serovar/biotype differentiation among major chicken isolates.

109  
110 Routine screening of flocks for the presence of *Salmonella* can be done by  
111 conventional serology which is expensive, as well as time- and labor-consuming.  
112 Based on the restricted number of major serovars found in chicken, extensive  
113 molecular techniques are not always cost-effective, and simpler more focused  
114 approaches could serve as rapid early diagnostic tests. Here, we took advantage of a  
115 small gap in gene *steB* of biotype Pullorum that was predicted by genomic analysis  
116 (18) to design primers that hybridize to biotype Gallinarum, but not Pullorum DNA,  
117 permitting a one-step PCR differentiation of the two biotypes (Table 2). This method  
118 shortened a previously described two-step technique (19). The addition of primers for  
119 additional chicken-associated serovars all in one multiplex PCR analysis is useful for  
120 the diagnosis of *Salmonella* in these birds. Although the designed probes are specific  
121 for the identification of serovars Heidelberg, Enteritidis and Gallinarum, serovar  
122 Kentucky shares its PCR profile with serovar Albany, which is not a major chicken  
123 isolate in the USA (13, 14). If needed, these two serovars could be differentiated by a  
124 flagellin-specific PCR (Suppl. Fig. 1). Finally, rarer serovars for which genomic data  
125 are currently unavailable might theoretically share one of the described PCR profiles,  
126 but as such serovars are significantly less frequent in chicken (13, 14), this would be

of minor concern.

Taken together, this study used (a) genomic sequence data for *Salmonella* to design a chicken-specific multiplex PCR diagnostic test and (b) an extensive library of *Salmonella* strains and serovars to validate the specificity of the method for the identification and differentiation of major avian-associated serovars. This simple and economical test should be useful for specific screening of poultry flocks, particularly for developing countries or backyard flocks and game birds in developed countries.

#### Acknowledgements

We thank Leon De Masi for critically reading the manuscript. This work was supported by China Scholarship Council grant [2013]3018 to CZ, and funds from NIH grant AI098041, USDA grant 2013-67015-21285 and the PennVet Center for Host-Microbial Interactions to DMS.

#### References

1. **Cohen ML, Tauxe RV.** 1986. Drug-resistant *Salmonella* in the United States: an epidemiologic perspective. *Science* **234**:964-969.
2. **Crim SM, Iwamoto M, Huang JY, Griffin PM, Gilliss D, Cronquist AB, Cartter M, Tobin-D'Angelo M, Blythe D, Smith K, Lathrop S, Zansky S, Cieslak PR, Dunn J, Holt KG, Lance S, Tauxe R, Henao OL, Centers for Disease C, Prevention.** 2014. Incidence and trends of infection with pathogens transmitted commonly through food--Foodborne Diseases Active Surveillance Network, 10 U.S. sites, 2006-2013. *MMWR Morb Mortal Wkly Rep* **63**:328-332.
3. **Foley SL, Nayak R, Hanning IB, Johnson TJ, Han J, Ricke SC.** 2011. Population dynamics of *Salmonella enterica* serotypes in commercial egg and poultry production. *Appl Environ Microbiol* **77**:4273-4279.
4. **Keller LH, Schifferli DM, Benson CE, Aslam S, Eckroade RJ.** 1997. Invasion of chicken reproductive tissues and forming eggs is not unique to *Salmonella enteritidis*. *Avian Dis* **41**:535-539.
5. **Keller LH, Benson CE, Krotec K, Eckroade RJ.** 1995. *Salmonella enteritidis* colonization of the reproductive tract and forming and freshly laid eggs of chickens. *Infect Immun* **63**:2443-2449.
6. **Foley SL, Johnson TJ, Ricke SC, Nayak R, Danzeisen J.** 2013. *Salmonella* pathogenicity and host adaptation in chicken-associated serovars. *Microbiol Mol Biol Rev* **77**:582-607.
7. **Shivaprasad HL, Methner U, Barrow PA.** 2013. *Salmonella* infections in the domestic fowl, p 162-192. In Barrow PA, Methner U (ed), *Salmonella* in domestic animals, 2nd ed. CABI, Wallingford, UK.
8. **Barrow PA, Freitas Neto OC.** 2011. Pullorum disease and fowl typhoid--new thoughts on old diseases: a review. *Avian Pathol* **40**:1-13.
9. **Andino A, Hanning I.** 2015. *Salmonella enterica*: survival, colonization, and

- 171 virulence differences among serovars. ScientificWorldJournal **2015**:520179.
- 172 10. **Baumler AJ, Hargis BM, Tsois RM.** 2000. Tracing the origins of *Salmonella*
- 173 outbreaks. Science **287**:50-52.
- 174 11. **Routh JA, Pringle J, Mohr M, Bidoi S, Arends K, Adams-Cameron M,**
- 175 **Hancock WT, Kissler B, Rickert R, Folster J, Tolar B, Bosch S, Barton**
- 176 **Behraves C, Williams IT, Gieraltowski L.** 2015. Nationwide outbreak of
- 177 multidrug-resistant *Salmonella* Heidelberg infections associated with ground
- 178 turkey: United States, 2011. Epidemiol Infect
- 179 doi:10.1017/S0950268815000497:1-8.
- 180 12. **Centers for Disease C, Prevention.** 2013. Outbreak of *Salmonella* Heidelberg
- 181 infections linked to a single poultry producer -- 13 states, 2012-2013. MMWR
- 182 Morb Mortal Wkly Rep **62**:553-556.
- 183 13. **Waltman D, Sellers R.** Report of the Committee on *Salmonella*, p. In (ed),
- 184 14. **Waltman D, Sellers R.** Report of the Committee on *Salmonella*, p. In (ed),
- 185 15. **Mulvey MR, Boyd DA, Finley R, Fakharuddin K, Langner S, Allen V, Ang L,**
- 186 **Bekal S, El Bailey S, Haldane D, Hoang L, Horsman G, Louis M, Robberts L,**
- 187 **Wylie J.** 2013. Ciprofloxacin-resistant *Salmonella enterica* serovar Kentucky in
- 188 Canada. Emerg Infect Dis **19**:999-1001.
- 189 16. **Rickert-Hartman R, Folster JP.** 2014. Ciprofloxacin-resistant *Salmonella*
- 190 *enterica* serotype Kentucky sequence type 198. Emerg Infect Dis **20**:910-911.
- 191 17. **Le Hello S, Hendriksen RS, Doublet B, Fisher I, Nielsen EM, Whichard JM,**
- 192 **Bouchrif B, Fashae K, Granier SA, Jourdan-Da Silva N, Cloeckert A,**
- 193 **Threlfall EJ, Angulo FJ, Aarestrup FM, Wain J, Weill FX.** 2011. International
- 194 spread of an epidemic population of *Salmonella enterica* serotype Kentucky
- 195 ST198 resistant to ciprofloxacin. J Infect Dis **204**:675-684.
- 196 18. **Yue M, Rankin SC, Blanchet RT, Nulton JD, Edwards RA, Schifferli DM.**
- 197 2012. Diversification of the *Salmonella* Fimbriae: A Model of Macro- and
- 198 Microevolution. PLoS ONE **7**:e38596.
- 199 19. **Cheraghchi N, Khaki P, Moradi Bidhendi S, Sabokbar A.** 2014. Identification
- 200 of Isolated *Salmonella enterica* Serotype gallinarum Biotype Pullorum and
- 201 Gallinarum by PCR-RFLP. Jundishapur J Microbiol **7**:e19135.
- 202
- 203

**Table 1.** List of primers and concentration used for PCR, with targeted DNA and amplicons sizes.

Primers	Sequences (5' to 3')	Final primer concentrations (pmole/ml)	Targeted genes or loci	Targeted DNA (species, serovars)	Amplicon sizes (bp)	Accession # and Nt segments
befC-F	GGGTGGCGGAAAACTATTTC	0.6	<i>befC</i>	<i>S. enterica</i>	993	AM933172 25665-26657
befC-R	CGGCACGCGGAATAGAGCAC					
heli-F	ACAGCCCCTGTTTAAATGGTG	2	orf (predicted helicase)	Heidelberg	782	CP005995 3226024-3226805
heli-R	CGCGTAATCGAGTAGTTGCC					
steB-F	TGTCGACTGGGACCCGCCGCCCGC	2	<i>steB</i>	Gallinarum biotype Gallinarum <sup>1</sup>	636	AM933173 2976016-2976651
steB-R	CCATCTTGTAGCGACCAT					
rhs-F	TCGTTTACGGCATTACACAAGTA	2.6	rhs locus	Gallinarum	402	AM933173 334109-334510
rhs-R	CAAACCCAGAGCCAATCTTATCT					
sdf-F	TGTGTTTATCTGATGCAAGAG	2.6	sdf locus	Enteritidis	293	AF370716 4950-5242
sdf-R	CGTTCTTCTGGTACTTCAGATGAC					
gly-F	TTCCAATTGAAACGAGTGCGG	2.6	orf "gly" (hypothetical protein)	Kentucky	170	ABE101000007 116981 - 117150
gly-R	ACTAACCCTGGGTGTGTGCTGT					

<sup>1</sup> Absent in biotype Pullorum (Accession number CP006575, locus\_tag="1137\_00945"), but also present in serovars Enteritidis, Heidelberg, Kentucky and group 1 serovars, as listed in Table 2.

211 **Table 2:** Bacterial strains used for confirm the specificity of the multiplex PCR assay.  
212

<i>Salmonella enterica</i> serovars and biotypes <sup>1</sup>	Multiplex PCR positive for					
	<i>bcfC</i>	<i>heli</i>	<i>steB</i>	<i>rhs</i>	<i>sdf</i>	<i>gly</i>
Heidelberg (2)	+	+	+	-	-	-
Enteritidis (11)	+	-	+	-	+	-
Kentucky (4)	+	-	+	-	-	+
Gallinarum biotype Gallinarum (16)	+	-	+	+	-	-
Gallinarum biotype Pullorum (7)	+	-	-	+	-	-
Others: Group 1(68) <sup>2</sup>	+	-	+	-	-	-
Others: Group 2 (20) <sup>3</sup>	+	-	-	-	-	-
Non <i>Salmonella</i> strains (5) <sup>4</sup>	-	-	-	-	-	-


213 <sup>1</sup> Numbers of strains in brackets (see Suppl. Table 1)

214 <sup>2</sup> Other *S. enterica* serovars (group 1) that have the same PCR profile: Paratyphi A (4 isolates),  
215 Paratyphi B var. Java (1), Agona (4), Abortusequi (2), Abortusovis (2), Saintpaul (3), Stanleyville  
216 (1), Typhisuis (2), Braenderup (5), Choleraesuis (24), Ohio (1), Thompson (1), Hadar (2),  
217 Muenchen (2), Newport (6), Berta (2), Dublin (2), Panama (1), Typhi (1), Agoueve (1), Cerro (1).

218 <sup>3</sup> Other *S. enterica* serovars (group 2) that have the same PCR profile: Schwarzengrund (3).  
219 Typhimurium (2), Bareilly (1), Hartford (1), Montevideo (2), Oranienburg (3), Javiana (6),  
220 Mississippi (1), Pomona (1).

221 <sup>4</sup> Three *E. coli* and two *Yersinia* strains.  
222  
223

224  
225  
226  
227  
228  
229  
230  
231  
232



233  
234  
235  
236  
237  
238  
239  
240  
241  
242  
243  
244  
245  
246  
247  
248

**Fig. 1.** Agarose gel (1.5%) of multiplex PCR amplicons from different bacterial strains. Representative gel from three comparable experiments. Lanes 1 and 10, 100 bp DNA ladder (NEB, Ipswich, MA); Lane 2, *Escherichia coli* (DH5a, negative control); Lane 3, *S. enterica* group 2, according to Table 2; Lane 4, *S. enterica* group 1, according to Table 2; Lane 5, *S. enterica* serovar Enteritidis; Lane 6, *S. enterica* serovar Heidelberg; Lane 7, *S. enterica* serovar Kentucky; Lane 8, *S. enterica* serovar Gallinarum biotype Pullorum; Lane 9, *S. enterica* serovar Gallinarum biotype Gallinarum.



