

# Resistance to carbapenems: mechanisms and detection methods

**Prof. Andrea Endimiani, MD, PhD**

Institute for Infectious Diseases (IFIK)  
University of Bern, Switzerland

1940

Penicillins

Broad-spectrum cephalosporins

Extended-spectrum cephalosporins

Carbapenems

Imipenem  
Meropenem  
Ertapenem  
Doripenem

2000

Penicillinases

Broad-spectrum  $\beta$ -lactamases

*S. aureus, H. influenzae, N. gonorrhoeae*  
*M. catarrhalis, V. cholerae, Enterobacteriaceae*  
*P. aeruginosa, Acinetobacter*

Extended-spectrum  $\beta$ -lactamases  
(ESBLs)

pAmpC  $\beta$ -lactamases

Carbapenemases


*Enterobacteriaceae*  
*P. aeruginosa*  
*Acinetobacter*

### *E. coli*, 2015: R to carbapenems




EU: 2012: 0.1% - 2015: 0.1%

### *K. pneumoniae*, 2015: R to carbapenems




EU: 2012: 6.2% - 2015: 8.1%

### *P. aeruginosa*, 2015: R to carbapenems



EU: 2012: 17.2% - 2015: 17.8%


### *Acinetobacter* spp., 2015: R to carbapenems



Italy: 2012: 83.1% - 2015: 78.3%

# How they become resistant to $\beta$ -lactams?


## Inactivation of $\beta$ -lactams




>2000 described  
(with very  $\neq$  spectrum)

Naturally present  
(e.g., *P. aeruginosa* are R ertapenem)

Acquired  
(chromosomal disruptions)



## Plasmid conjugation



Exchange and spread among  
different bacterial species


Naturally present  
(e.g., *P. aeruginosa*)

Acquired  
(chromosomal overexpression)


Development of ertapenem resistance in a patient with mediastinitis caused by *Klebsiella pneumoniae* producing an extended-spectrum  $\beta$ -lactamase

David Skurnik,<sup>1†</sup> Sigismond Lasocki,<sup>2</sup> Sylvie Bremont,<sup>3</sup>  
 Claudette Muller-Serieys,<sup>1</sup> Marie Dominique Kitzis,<sup>4</sup> Patrice Courvalin,<sup>3</sup>  
 Antoine Andremon<sup>1</sup> and Philippe Montravers<sup>2</sup> JMM, 2010

Triple coronary  
bypass



ERTA: 1 g/day  
 IMP: 4 g/day



Identical by PFGE

IMP: from  $\leq 0.125$  to 8 mg/L  
 MEM: from  $\leq 0.125$  to 12 mg/L  
 ERTA: from 0.094 to 8 mg/L

# MINIREVIEW

ANTIMICROBIAL AGENTS AND CHEMOTHERAPY, Mar. 2010, p. 969–976

## Updated Functional Classification of $\beta$ -Lactamases<sup>v</sup>

Karen Bush<sup>1\*</sup> and George A. Jacoby<sup>2</sup>

TABLE 1. Classification schemes for bacterial  $\beta$ -lactamases, expanded from Bush et al. (16)

Bush-Jacoby group (2009)	Bush-Jacoby-Medeiros group (1995)	Molecular class (subclass)	Distinctive substrate(s)	Inhibited by		Defining characteristic(s)	Representative enzyme(s)
				CA or TZB <sup>a</sup>	EDTA		
1	1	C	Cephalosporins	No	No	Greater hydrolysis of cephalosporins than benzylpenicillin; hydrolyzes cephamycins	<i>E. coli</i> AmpC, P99, ACT-1, CMY-2, FOX-1, MIR-1
1e	NI <sup>b</sup>	C	Cephalosporins	No	No	Increased hydrolysis of ceftazidime and often other oxymino- $\beta$ -lactams	GC1, CMY-37
2a	2a	A	Penicillins	Yes	No	Greater hydrolysis of benzylpenicillin than cephalosporins	PC1
2b	2b	A	Penicillins, early cephalosporins	Yes	No	Similar hydrolysis of benzylpenicillin and cephalosporins	TEM-1, TEM-2, SHV-1
2be	2be	A	Extended-spectrum cephalosporins, monobactams	Yes	No	Increased hydrolysis of oxymino- $\beta$ -lactams (cefotaxime, ceftazidime, ceftriaxone, cefepime, aztreonam)	TEM-3, SHV-2, CTX-M-15, PER-1, VEB-1
2br	2br	A	Penicillins	No	No	Resistance to clavulanic acid, sulbactam, and tazobactam	TEM-30, SHV-10
2ber	NI	A	Extended-spectrum cephalosporins, monobactams	No	No	Increased hydrolysis of oxymino- $\beta$ -lactams combined with resistance to clavulanic acid, sulbactam, and tazobactam	TEM-50
2c	2c	A	Carbenicillin	Yes	No	Increased hydrolysis of carbenicillin	PSE-1, CARB-3
2ce	NI	A	Carbenicillin, cefepime	Yes	No	Increased hydrolysis of carbenicillin, cefepime, and ceftiofime	RTG-4
2d	2d	D	Cloxacillin	Variable	No	Increased hydrolysis of cloxacillin or oxacillin	OXA-1, OXA-10
2de	NI	D	Extended-spectrum cephalosporins	Variable	No	Hydrolyzes cloxacillin or oxacillin and oxymino- $\beta$ -lactams	OXA-11, OXA-15
2df	NI	D	Carbapenems	Variable	No	Hydrolyzes cloxacillin or oxacillin and carbapenems	OXA-23, OXA-48
2e	2e	A	Extended-spectrum cephalosporins	Yes	No	Hydrolyzes cephalosporins. Inhibited by clavulanic acid but not aztreonam	CepA
2f	2f	A	Carbapenems	Variable	No	Increased hydrolysis of carbapenems, oxymino- $\beta$ -lactams, cephamycins	KPC-2, IMI-1, SME-1
3a	3	B (B1)	Carbapenems	No	Yes	Broad-spectrum hydrolysis including carbapenems but not monobactams	IMP-1, VIM-1, CcrA, IND-1
		B (B3)					
3b	3	B (B2)	Carbapenems	No	Yes	Preferential hydrolysis of carbapenems	L1, CAU-1, GOB-1, FEZ-1, CphA, Sfh-1
NI	4	Unknown					

<sup>a</sup> CA, clavulanic acid; TZB, tazobactam.

<sup>b</sup> NI, not included.

OXA-24/40

NDMs



Occurrence of carbapenemase-producing *Klebsiella pneumoniae* and *Escherichia coli* in the European survey of carbapenemase-producing Enterobacteriaceae (EuSCAPE): a prospective, multinational study Grundmann et al., LID, Feb 2017

Nov, 2013 - Apr 2014  
455 hospitals  
36 countries

<i>Klebsiella pneumoniae</i>			<i>Escherichia coli</i>			Sentinel hospitals (mean beds†)	Incidence per 10000 admissions‡		Incidence per 100000 patient-days§	
Submitted non-susceptible isolates (n)	Confirmed carbapenemase-producing isolates* (n)	Comparator isolates (n)	Submitted non-susceptible isolates (n)	Confirmed carbapenemase-producing isolates* (n)	Comparator isolates (n)		Rate (hospitals)	Rank	Rate (hospitals)	Rank
1203	850 70.6%	1098	194	77 39.7%	208	455 (800)	1.3 (321)	..	2.51 (268)	..

Hospitals submitting carbapenem non-susceptible <i>E coli</i> isolates (n)	Number of submitted carbapenem non-susceptible <i>E coli</i> isolates (n)	Confirmed carbapenemase-producing <i>E coli</i> isolates					Other (n, %)*
		KPC (n, %)	NDM (n, %)	OXA-48-like	VIM (n, %)	Total (n, %)	
105	194	14 (7.2)	20 (10.3)	43 (22.2)	0	77 (39.7)	117 (60.3)
		18.2%	26.0%	55.8%			

Hospitals submitting carbapenem non-susceptible <i>K pneumoniae</i> isolates (n)	Number of submitted carbapenem non-susceptible <i>K pneumoniae</i> isolates	Confirmed carbapenemase-producing <i>K pneumoniae</i> isolates					Other (n, %)*
		KPC (n, %)	NDM (n, %)	OXA-48-like (n, %)	VIM (n, %)	Total (n, %)	
251	1203	379 (31.5)	93 (7.7)	310 (25.8)	68 (5.7)	850 (70.7)	353 (29.3)
		44.6%	10.9%	36.5%	0.08%		

1. KPC: 42.4%
2. OXA-48-like: 38.1%
3. NDM: 12.2%
4. VIM: 7.3%

Good hydrolysis → "No" hydrolysis

Penicillins

NS-Ceph

3GCs

4GCs

Anti-MRSA

Monobactams

Cefamycins

Carbapenems

Imipenem  
Meropenem  
Ertapenem  
Doripenem



**ESBLs**  
(Class A; e.g., CTX-Ms)



**KPCs**  
(Class A; e.g., KPC-2, KPC-3)



**Metallo-β-lactamases, MBLs**  
(Class B; e.g., NDM-1, VIM, IMP)



**AmpCs**  
(Class C; e.g., CMYs)



**OXAs (broad/ext.-spectrum)**  
(Class D)



**OXA carbapenemases**  
(Class D; e.g., OXA-23/-24/-48)

\*Catalytic efficiency for imipenem:  $k_{cat}/K_m$  [ $\mu\text{M}^{-1}/\text{s}^{-1}$ ]




World map showing the global distribution of KPC-2 to KPC-22. The map is color-coded: red for endemic spread, dark blue for sporadic spread, light blue for recorded, and white for not recorded. Numbers are placed within various countries to indicate specific data points.

Legend:


- Endemic spread of KPC producers
- Sporadic spread of KPC producers
- KPC recorded
- Not recorded

**Tn4401 is carried by MDR plasmids (IncF, I2, X, A/C, R, ColE1) of which some hyperepidemic**



MLST	<i>gapA</i>	<i>infB</i>	<i>mdh</i>	<i>pgi</i>	<i>phoE</i>	<i>rpoB</i>	<i>tonB</i>
ST11	3	3	1	1	1	1	4
ST37	2	9	2	1	13	1	16
ST258	3	3	1	1	1	1	79
ST277	3	1	1	1	1	1	43
ST307	4	1	2	52	1	1	7
ST327	2	1	1	1	10	1	19
ST340	3	3	1	1	1	1	18
ST376	2	6	1	3	8	1	27
ST384	18	23	56	63	80	43	51
ST388	16	24	59	27	29	22	105
ST512	54	3	1	1	1	1	79


# OXA-48-like carbapenemases (class D)



## OXA-48-like carbapenemases: the phantom menace

*J Antimicrob Chemother* 2012; **67**: 1597–1606

Laurent Poirel\*, Anaïs Potron and Patrice Nordmann



OXA-48 *K. pneumoniae*


ESBL/OXA-48 *K. pneumoniae*

OXA-48 variants: OXA-162, -163, -181, -204, -232, -244, -245

## Common hyperepidemic L plasmid (non-MDR)

Strain	Country	Size (kb) <sup>a</sup>	Plasmid feature				Plasmid-associated resistance marker
			RepA	TraU	ParA	Inc group	
<i>K. pneumoniae</i> 11978	Turkey	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> Bel	Belgium	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> Lib	Lebanon	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> 8	Turkey	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> Egy	Egypt	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> 5A	Turkey	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> 7A	Turkey	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> Bey	Lebanon	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> 3A	Turkey	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> 4A	Turkey	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> 17A	Turkey	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> Bou	France	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> Dia	France	ca. 62	+	+	+	L/M	None
<i>K. pneumoniae</i> C2	Morocco	ca. 62	+	+	+	L/M	None
<i>K. oxytoca</i> A7	Morocco	ca. 62	+	+	+	L/M	None
<i>E. coli</i> 1	Turkey	ca. 62	+	+	+	L/M	None
<i>E. cloacae</i> 501	France	ca. 62	+	+	+	L/M	None
<i>E. cloacae</i> Bou	Morocco	ca. 62	+	+	+	L/M	None
<i>E. cloacae</i> D4	Morocco	ca. 62	+	+	+	L/M	None

## Composite transposon Tn1999




## Hyper epidemic clones

*E. coli* ST131

*K. pneumoniae* ST11, ST147

# NDM carbapenemases (class B, MBL)



Characterization of a New Metallo- $\beta$ -Lactamase Gene, *bla*<sub>NDM-1</sub>, and a Novel Erythromycin Esterase Gene Carried on a Unique Genetic Structure in *Klebsiella pneumoniae* Sequence Type 14 from India<sup>7</sup>

Donggun Yong,<sup>1,2</sup> Mark A. Toleman,<sup>2</sup> Christian G. Giske,<sup>3</sup> Hyun S. Cho,<sup>4</sup> Kristina Sundman,<sup>5</sup> Kyungwon Lee,<sup>1</sup> and Timothy R. Walsh<sup>2\*</sup>

ANTIMICROBIAL AGENTS AND CHEMOTHERAPY, Dec. 2009, p. 5046–5054

2008: male with UTI/colonized (back from India)

Emergence of a new antibiotic resistance mechanism in India, Pakistan, and the UK: a molecular, biological, and epidemiological study

*Lancet Infect Dis* 2010;

2009: many UK patients who traveled in India/Pakistan

## Very promiscuous *bla* gene(s)

### Many species

2009–2012 (950 NDM strains)

<i>K. pneumoniae</i>	37.8%
<i>E. coli</i>	28.2%
Other Gram-neg. bacilli	12.0%
Other <i>Enterobacteriaceae</i>	11.8%
<i>E. cloacae</i>	5.5%
<i>A. baumannii</i>	3.8%
<i>P. aeruginosa</i>	1.0%


Not necessarily hyperepidemic clones!

### Many MDR plasmids


IncA/C,  
IncF, IncR,  
IncH, IncN,  
IncX, IncL/M

Frequently co-carrying important ARGs (e.g., 23S rRNA methylases)


### Genetic structures



# Standard culture approach



Even 5–6 days  
to get the final report



Empirical  
treatment

Isolation



Directed  
treatment

## Septic shock

*Inappropriate antibiotics within the first 6 hours is associated with a 5-fold higher mortality risk*

[Kumar A *et al.*, Chest, 2009]

The Changing Role of the  
Clinical Microbiology  
Laboratory in Defining Resistance  
in Gram-negatives IDCNA, 2016

Andrea Endimiani, MD, PhD<sup>a,\*</sup>, Michael R. Jacobs, MD, PhD<sup>b</sup>


# CLSI 2017: M100-S27

Laboratories using *Enterobacteriaceae* MIC breakpoints for carbapenems described in M100-S20 (January 2010) should perform the MHT, the Carba NP test, mCIM, and/or a molecular assay when isolates of *Enterobacteriaceae* are suspicious for carbapenemase production based on imipenem or meropenem MICs of 2–4 µg/mL or ertapenem MIC of 2 µg/mL (refer to Tables 3B, 3C, and 3D). After implementation of the current breakpoints, these additional tests do not need to be performed other than for epidemiological or infection control purposes (refer to Table 3B).

	CLSI, Jan 2010		CLSI, 2017		EUCAST, 2017	
	S	R	S	R	S [screening]	R
IMP	≤4	≥16	≤1	≥4	≤2 [≥1]	≥16
MEM	≤4	≥16	≤1	≥4	≤2 [≥0.12]	≥16

## Screening tests

MHT (Table 3B)	Carba NP (Table 3C)	mCIM (Table 3D)
<i>Enterobacteriaceae</i> that are not susceptible to one or more carbapenems	<i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , and <i>Acinetobacter</i> spp. that are not susceptible to one or more carbapenems	<i>Enterobacteriaceae</i> that are not susceptible to one or more carbapenems



What Is the Appropriate Meropenem MIC for Screening of Carbapenemase-Producing *Enterobacteriaceae* in Low-Prevalence Settings?  
Fattouh R et al., AAC, March 2016

1,022 *Enterobacteriaceae* (189 CPE)

### CLSI:

- 14% had MICs in the S range
- Screening tests (if used) captured 86% of the CPE

### EUCAST:

- 20% had MICs in the S range
- Screening tests captured 98% of the CPE




# CLSI: Modified Hodge test


**Cloverleaf test (modified Hodge test) for detecting carbapenemase production in *Klebsiella pneumoniae*: be aware of false positive results**

*J Antimicrob Chemother* 2010; 65: 249–251

Cecilia G. Carvalhaes\*, Renata C. Picão, Adriana G. Nicoletti, Danilo E. Xavier and Ana C. Gales



This test indicates the possible production of enzyme(s) able to hydrolyze carbapenems (i.e., meropenem) = generic CARBAPENEMASE



**Figure 1.** Phenotypic carbapenemase detection by MHT applying the inoculum recommended by CLSI for ertapenem (a) and meropenem discs (b). 1, *K. pneumoniae* ATCC BAA-1705, positive result; 2, *K. pneumoniae* ATCC BAA-1706, negative result; 3, CTX-M-producing *K. pneumoniae* clinical isolate; 4, KPC-producing *K. pneumoniae* clinical isolate. Arrows indicate the similar size of *E. coli* ATCC 25922 grown within the carbapenem disc inhibition zones when testing both carbapenemase producer and carbapenemase non-producer isolates.

ESBL or pAmpC producers  
with porin(s) loss

Not good for NDMs


# CLSI: Modified Carbapenem Inactivation Method (mCIM)



1  $\mu$ L loop bacteria  
2 mL TSB  
10  $\mu$ g MEM




4 h incubation



ATCC 25922 *E. coli*  
0.5 McFarland  
MH agar



35°C $\pm$ 2  
18-24 h



6-15 mm: Carbapenemase positive  
16-18 mm: Indeterminate results  
 $\geq$  19 mm: No carbapenemases

**Modified Carbapenem Inactivation  
Method for Phenotypic Detection of  
Carbapenemase Production among  
*Enterobacteriaceae*** Pierce VM, JCM, Aug 2017

117 *Enterobacteriaceae* (94 CPE)


Single clinical laboratory  
**Sens: 99%; Spec: 100%**

Nine laboratory study  
**Sens: 97%; Spec: 99%**



# Rapid Detection of Carbapenemase-producing *Enterobacteriaceae*

Nordmann, EID, 2012



CLSI: NP test

<2 hrs

Most CPE turned in 10-30 min

!!!! OXA-48-like  
not always detected

Need pre-extraction of enzymes (~30 min)

Blue-Carba, an Easy Biochemical Test for Detection of Diverse Carbapenemase Producers Directly from Bacterial Cultures  
J. Pires, A. Novals, L. Peixe December 2013 Volume 51 Number 12 Journal of Clinical Microbiology


Bromothymol blue

A: NDM-1 *E. coli*

B: OXA-23 *A. baumannii*

C: OXA-48 *K. pneumoniae*

D: *E. coli* ATCC 25922



Direct use of colonies!

# Non-molecular detection of carbapenemases in Enterobacteriaceae clinical isolates

A. Aguirre-Quinonero<sup>a</sup>, L. Martínez-Martínez<sup>a, b, \*</sup> J Infect Chemother 23 (2017) 1–11


## Inhibitor based approaches

	Synergy with MPM			
	DPA EDTA	PBA	DPA + PBA	Cloxacillin
KPC	–	+	–	–
MBL	+	–	–	–
OXA-48-like	–	–	–	–
KPC + MBL	v	v	+	–
AmpC + porin loss	–	+	–	+
ESBL + porin loss	–	–	–	–


MPM: Meropenem, DPA: dipicolinic acid, PBA: phenyl boronic acid.

## Combination disk test (CDT)

Rosco Diagnostica KPC/MBL confirmation kit




KPC-producing *K. pneumoniae*



NDM-producing *K. pneumoniae*

## Gradient diffusion strips

Etest: Imipenem+EDTA / Imipenem



VIM-producing *K. pneumoniae*

# Evaluation of a DNA microarray for rapid detection of the most prevalent extended-spectrum $\beta$ -lactamases, plasmid-mediated cephalosporinases and carbapenemases in Enterobacteriaceae, *Pseudomonas* and *Acinetobacter*

Bogaerts P, IJAA, 2016



50-85 Euro


~6 hrs

**ESBLs:** TEM, SHV, CTX-M, BEL, PER, GES, VEB

**pAmpCs:** CMY, DHA, FOX, ACC-1, ACT/MIR

**Carbapenemases:** KPC, NDM, VIM, IMP, GIM, SPM, OXA-48, OXA-23, OXA-24, OXA-58

HybC	TEM 104E (1)	TEM 104K (2)	DNA-C (3)	TEM 164R (4)	TEM 164S (5)		TEM 164H (7)	IMP (8)	IMP (9)
CTR (10)	TEM 238G (11)	TEM 238S (12)	DNA-C (13)	SHV-all (14)	SHV 238G (15)	SHV G238S (16)	SHV G238A (17)	SHV 240E (18)	SHV E240K (19)
CTR (20)	OXA-48 (21)	OXA-48 (22)	negC (23)	VIM (24)	VIM (25)	CTX-M1 fam (26)	CTX-M2 fam (27)	CTX-M9 fam (28)	CTX-M8+25 fam (29)
VIM-all (30)	KPC (31)	KPC (32)	HybC	FOX (33)	CTX-M1-A (34)	CTX-M1-B (35)	VIM-all (36)	CTX-M1-C (37)	CTX-M1-D (38)
PER I (39)	NDM (40)	PER I (41)	PER II (42)	TEM -all (43)	TEM 164S (44)	TEM 164C (45)	TEM 164H (46)	CMY II (47)	CMY II (48)
HybC	NDM (49)	VEB (50)	VEB (51)	SHV-all (52)	SHV 238G (53)	SHV G238S (54)	SHV G238A (55)	SHV 240E (56)	SHV E240K (57)
GIM (58)	DHA (59)	ACT/MIR (60)	GIM (61)	CMY I / MOX (62)	ACC (63)	SPM (64)	SPM (65)	BEL (66)	BEL (67)
OXA-23 (68)	OXA-23 (69)	TEM 104K (70)	DNA-C (71)	OXA-24 (72)	OXA-24 (73)	OXA-58 (74)	OXA-58 (75)	IMP (76)	IMP (77)
negC (78)	GES-all (79)	TEM 238S (80)	HybC	GES CARBA (81)	GES ESBL (82)	<i>mcr</i>	<i>mcr</i>	<i>mcr</i>	
	DHA (88)	ACT/MIR (89)	DNA-C (90)	CMY I / MOX (91)	ACC (92)	CTX-M1 fam (93)	CTX-M2 fam (94)	CTX-M9 fam (95)	CTX-M8+25 fam (96)



Sensitivity and Specificity: both ~100%

# Amplex Eazyplex




- Preparation (**5 min**)
- No DNA extraction
- Run time (**<30 min**)

~55 Euro

## LAMP

Loop-mediated isothermal Amplification  
Real-time fluorescent measurement



## eazyplex® SuperBug CRE

KPC,  
NDM,  
OXA-48  
and OXA 181,  
VIM, as well as  
CTX-M-1 and  
CTX-M-9 group

## eazyplex® SuperBug mcr-1

## eazyplex® SuperBug complete

### Version A

KPC  
NDM  
OXA-48  
VIM  
OXA-23 group  
OXA-40 group  
**OXA-58 group**

### Version B

KPC  
NDM  
OXA-48  
VIM  
OXA-23 group  
OXA-40 group  
**OXA-181**

## Evaluation of the eazyplex® SuperBug CRE system for rapid detection of carbapenemases and ESBLs in clinical Enterobacteriaceae isolates recovered at two Spanish hospitals

J Antimicrob Chemother 2015; 70: 1047 – 1050

94 carbapenemase-producing *Enterobacteriaceae*



**100% agreement** with PCR/Sequencing of *bla* genes

KPC-2	<i>K. pneumoniae</i> (2), <i>Enterobacter cloacae</i> (2), <i>E. coli</i> (1)
KPC-3	<i>K. pneumoniae</i> (22), <i>E. cloacae</i> (1), <i>E. coli</i> (2)
KPC-3 + CTX-M-15	<i>K. pneumoniae</i> (1)
VIM-1	<i>K. pneumoniae</i> (2), <i>Klebsiella oxytoca</i> (1), <i>E. coli</i> (2), <i>E. cloacae</i> (3), <i>Citrobacter freundii</i> (1), <i>Serratia marcescens</i> (2), <i>Raoultella ornithinolytica</i> (1)
VIM-1 + CTX-M-10	<i>K. oxytoca</i> (1)
VIM-1 + CTX-M-15	<i>K. pneumoniae</i> (1), <i>K. oxytoca</i> (1)
VIM-1 + CTX-M-32	<i>E. coli</i> (1)
VIM-1 + CTX-M-14	<i>K. oxytoca</i> (1)
NDM-1	<i>Providencia rettgeri</i> (1)
NDM-5	<i>E. coli</i> (1)
OXA-48	<i>K. pneumoniae</i> (3), <i>K. oxytoca</i> (1), <i>E. coli</i> (2), <i>E. cloacae</i> (1), <i>Citrobacter koseri</i> (1), <i>Citrobacter braakii</i> (1)
OXA-48 + CTX-M-15	<i>K. pneumoniae</i> (30), <i>E. coli</i> (2), <i>Enterobacter aerogenes</i> (1), <i>Citrobacter amalonaticus</i> (1)
OXA-48 + CTX-M-14	<i>E. coli</i> (1)

# Rapid detection and identification of strains carrying carbapenemases directly from positive blood cultures using MALDI-TOF MS

Hoyos-Mallecot *et al.*, J Micro Meth; 105, 2014

19 carba producers (21 controls)  
Sens and Spec: 100% and 90%



0.30 Euro\*



Generic activity in ~4.5 hrs  
Can be reduced ~30 min-2 hrs



# Rapid detection of carbapenemase-producing Enterobacteriaceae from blood cultures

L. Dortet<sup>1</sup>, L. Brécard<sup>1</sup>, L. Poirel<sup>1,2</sup> and P. Nordmann<sup>1,2</sup>

*Clin Microbiol Infect* 2014; **20**: 340–344



The RAPIDEC Carba NP <sup>®</sup>Test

Carbapenemase types	Tested isolates (n)	Carba NP test on positive blood culture					
		Positive		Negative		Sensitivity (%)	Specificity (%)
		n	%	n	%		
KPC	50	50	100	0	0	100	100
IMP	27	27	100	0	0	100	100
VIM	37	37	100	0	0	100	100
NDM	33	33	100	0	0	100	100
OXA-48-like	46	42	91.3	4	8.7	91.3	100
No carbapenemase	74	0	0	74	100	97.9	—
Total results							100

Preparation from +BCs: ~3 hrs

Reaction in 5 min-1 h (max 2 hrs)


Sens and Spec 98% and 100%

# Check-Direct Screening, Check-Points

**Rectal  
swab**  
(Copan, ESwab)




- Real-time multiplex PCR
- Rapid preparation
- Time to results (<3 hrs)



**BD MAX**



**Reagent strip**



Check-Direct **CPE**  
Screen for **BD MAX™**

~30 Euro

KPC  
OXA-48-like  
VIM  
NDM



Performance of the BD MAX™ instrument with Check-Direct CPE real-time PCR for the detection of carbapenemase genes from rectal swabs, in a setting with endemic dissemination of carbapenemase-producing *Enterobacteriaceae* Antonelli *et al.*, DMID, 2016

Florence, Italy  
557 rectal swabs

ChromID CARBA SMART (w/wo broth enrichment)  
[5 samples not detected]

### Summary of the results

		Colonization status	
		Positive (n. 29)	Negative (n. 528)
CHSM and BERM	Positive	24 (82.8 %) <sup>a</sup>	0
	Negative	5 (17.2 %)	528 (100 %)
CPMix	Positive	29 (100 %)	15 (2.9 %)
	Negative	0	507 (96.0 %)
	Unresolved	-	6 (1.1 %)

### Limit of detection (LOD)

Target	LOD CPMix (CFU/ml)
<i>bla</i> <sub>KPC</sub>	9×10 <sup>3</sup>
<i>bla</i> <sub>VIM</sub>	4.5×10 <sup>2</sup>
<i>bla</i> <sub>OXA-48</sub>	8.5×10 <sup>1</sup>
<i>bla</i> <sub>NDM</sub>	7.3×10 <sup>1</sup>

23 KPC  
5 VIM  
1 OXA-48

TAT from 18-24 h (direct culture) or 48 h (broth enrichment)  
to 3 hrs

# Cepheid GeneXpert



- Add aliquot to elution, vortex, transfer to port 5
- Insert cartridge to station (overall, **1 min**)
- Run time (**<1 h**)

## Real-time multiplex PCR

- Smart fluidic system
- Filtering and Sonication (DNA)
- Fluorescent-labeled hybr. probes (6 colors)
- Internal control

**Xpert® Carba-R**

KPC  
NDM  
OXA-48-like  
VIM  
IMP-1

~50 Euro

# Multisite Evaluation of Cepheid Xpert Carba-R Assay for Detection of Carbapenemase-Producing Organisms in Rectal Swabs JCM, 54:7; 2016

M. Tato,<sup>a</sup> P. Ruiz-Garbajosa,<sup>a</sup> M. Traczewski,<sup>b</sup> A. Dodgson,<sup>c</sup> A. McEwan,<sup>c</sup> R. Humphries,<sup>d</sup> J. Hindler,<sup>d</sup> J. Veltman,<sup>e</sup> H. Wang,<sup>f</sup> R. Cantón<sup>a</sup>

4 centers (2 USA, 1 UK, 1 Spain)

July 2013 - Feb 2014

633 samples

50% of them spiked at  $1.1 \times 10^2$  to  $1.2 \times 10^3$  CFU/swab (LOD of Xpert)

## Results (obtained in 32-48 min)

Xpert Carba-R assay result	Clinical specimens (n = 383)	Contrived specimens (n = 250)	All specimens (n = 633)
Positive (single and/or combined targets)	42	107	149 23.5%
IMP-1	0	25	25
VIM	2	24	26
NDM	2	23	25
KPC	13	19	32
OXA-48	20	15	35
VIM + OXA-48	4	0	4
NDM + KPC	1	0	1
IMP-1 + NDM	0	1	1
Negative	341	143	484

## Performance for different targets


Target gene	Sensitivity (% [95% CI])	Specificity (% [95% CI])	PPV (%)	NPV (%)
IMP-1	96.3 (81.0–99.9)	100 (99.4–100)	100	99.8
VIM	93.5 (78.6–99.2)	99.8 (99.1–100)	96.7	99.7
NDM	100 (86.8–100)	99.8 (99.1–100)	96.3	100
KPC	96.7 (82.8–99.9)	99.3 (98.3–99.8)	87.9	99.8
OXA-48	95.0 (83.1–99.4)	99.8 (99.1–100)	97.4	99.7

Sensitivity: 96.6% - Specificity: 98.6%

# Verigene, Nanosphere



- Load cartridge, consumables, and sample (5 min)
- Automated sample preparation and processing
- Place slide from cartridge in reader (2.5 hrs)



**Microarray** approach by using Au-nanoprobe as reporter and silver reduction to enhance signal

## Gram-negatives cartridge

Species	Genus	Resistance
<i>Escherichia coli</i> *	<i>Acinetobacter</i> spp.	CTX-M (ESBL)
<i>Klebsiella pneumoniae</i>	<i>Citrobacter</i> spp.	IMP (carbapenemase)
<i>Klebsiella oxytoca</i>	<i>Enterobacter</i> spp.	KPC (carbapenemase)
<i>Pseudomonas aeruginosa</i>	<i>Proteus</i> spp.	NDM (carbapenemase)
<i>Serratia marcescens</i>		OXA (carbapenemase)
		VIM (carbapenemase)

~50 Euro

# Performance Evaluation of the Verigene Gram-Positive and Gram-Negative Blood Culture Test for Direct Identification of Bacteria and Their Resistance Determinants from Positive Blood Cultures in Hong Kong

Gilman K.H. Siu *et al.*, Plos One, Oct 2015

Multicenter study (4 hospitals)  
364 BCs (114 Gram-pos; 250 Gram-neg)  
Jan 2014 - May 2014

## ID for Gram-negatives (agreement 90.5%)

Organisms	No. (%) of isolates					Sensitivity (%)	Specificity (%)
	Total	Correctly identified	Not detected	Misidentified	No Call		
<i>E.coli</i>	165 (59.1)	158 (95.8)	6 (3.6)	1 <sup>c</sup>	-	95.8	100
<i>K. pneumoniae</i>	52 (18.6)	36 (69.2)	14 (27) <sup>d</sup>	2 <sup>e</sup>	-	69.2 <sup>d</sup>	100
<i>P. aeruginosa</i>	17 (6.1)	13 (76.5)	3 (17.6)	-	1	81.3	100
<i>Proteus spp.</i>	10 (3.6)	9 (90)	-	-	1	100	100
<i>Enterobacter spp.</i>	10 (3.6)	8 (80)	-	-	2	100	99.3
<i>Acinetobacter spp.</i>	4 (1.4)	4 (100)	-	-	-	100	100
<i>K. oxytoca</i>	3 (1.1)	-	3 (100)	-	-	0	99.3

*K. variicola* 80%

## 40-99 h faster than routine

Organisms	No. of Isolates <sup>b</sup>	Δ Time to Result <sup>a</sup>		
		Average Time to identification by Culture-Based Method (h)	Average Time to Result by Verigene Test (h)	Average (h)
<i>Staphylococcus spp.</i>	23	63.76	2.35	61.41
MSSA	7	49	2.35	46.65
MRSA	10	93.45	2.35	91.1
CNS	8	57.19	2.35	54.84
<i>Streptococcus spp.</i>	13	63.86	2.35	61.51
<i>S. pneumoniae</i>	3	42.83	2.35	40.48
β-haemolytic Strept.	4	71.62	2.35	69.27
Viridians Group	6	55.08	2.35	52.73
<i>Enterococcus spp.</i>	2	78.5	2.35	76.15
VSE	1	55.5	2.35	53.15
VRE	1	101.5	2.35	99.15
Enterobacteriaceae	79	45.98	1.88	44.1
<i>P. aeruginosa</i>	7	52.51	1.88	50.63
<i>Acinetobacter spp.</i>	1	76	1.88	74.12



## Drug-Resistant Organisms

Drug resistant Organisms	No. of isolates				Sensitivity	Specificity (%)
	Total	Correctly Detected	Not Detected	No Call		
Gram Positive						
MRSA	27	26	0	1	100	100
MRSE	4	4	0	0	100	100
VRE	1	1	0	0	100	100
Gram Negative						
Cefotaxime resistant <i>Enterobacteriaceae</i> (including ESBL producers)	61	38	22	1	63.3	100
ESBL producing <i>Enterobacteriaceae</i>	46	38	7	1	84.4	100
MDR <i>Acinetobacter</i>	3	3	0	0	100	100
Carbapenem-resistant <i>Pseudomonas</i>	2	0	2	0	0	100
Total	98	72	24	2	75.0	100


# BioFire FilmArray



- Preparation of the pouch
- Add pouch to FilmArray station (overall, **2 min**)
- Run time of about **1 h**



Magnetic beads



Multiplex PCR  
(RT-PCR for RNA target)

Multiplex real-time PCR  
and Melting analysis

FilmArray® BCID Panel targets:

~100 Euro

Gram+ Bacteria	Gram- Bacteria
<i>Enterococcus</i> <i>Listeria monocytogenes</i> <i>Staphylococcus</i> <i>Staphylococcus aureus</i> <i>Streptococcus</i> <i>Streptococcus agalactiae</i> <i>Streptococcus pyogenes</i> <i>Streptococcus pneumoniae</i>	<i>Acinetobacter baumannii</i> <i>Haemophilus influenzae</i> <i>Neisseria meningitidis</i> <i>Pseudomonas aeruginosa</i> <i>Enterobacteriaceae</i> <i>Enterobacter cloacae</i> complex <i>Escherichia coli</i> <i>Klebsiella oxytoca</i> <i>Klebsiella pneumoniae</i> <i>Proteus</i> <i>Serratia marcescens</i>
Yeast	Antibiotic Resistance
<i>Candida albicans</i> <i>Candida glabrata</i> <i>Candida krusei</i> <i>Candida parapsilosis</i> <i>Candida tropicalis</i>	<i>mecA</i> - methicillin resistance <i>vanA/B</i> - vancomycin resistance <i>KPC</i> - carbapenem resistance



## Evaluation of the FilmArray Blood Culture Identification Panel: Results of a Multicenter Controlled Trial

Salimnia H.. *et al.*, J Clin Microbiol; 54:3, 2016

8 centers in USA  
July 2012 - Feb 2014  
2,207 BC samples  
(processed within 8 h)

Species ID	Isolates detected <sup>a</sup> : BCID/comparator		No. of results: BCID/comparator				Sensitivity or PPA <sup>b</sup> : TP/(TP + FN) (%)	Specificity or NPA <sup>b</sup> : TN/(TN + FP) (%)
	Clinical arm	Seeded arm	TP +/+	FP +/-	FN -/+	TN -/-		
Gram-positive bacteria								
<i>Enterococcus</i>	102/101	29/29	127	4	3	2,073	127/130 (97.7)	2,073/2,077 (99.8)
<i>L. monocytogenes</i>	0/0	36/36	36	0	0	2,171	36/36 (100)	2,171/2,171 (100)
<i>Staphylococcus</i>	780/797	2/1	770	12	28	1,397	770/798 (96.5)	1,397/1,409 (99.1)
<i>S. aureus</i>	257/257	0/0	253	4	4	1,946	253/257 (98.4)	1,946/1,950 (99.8)
<i>Streptococcus</i>	140/141	63/62	198	5	5	1,999	198/203 (97.5)	1,999/2,004 (99.8)
<i>S. agalactiae</i> (group B)	18/18	18/18	36	0	0	2,171	36/36 (100)	2,171/2,171 (100)
<i>S. pneumoniae</i>	26/25	12/12	36	2	1	2,168	36/37 (97.3)	2,168/2,170 (99.9)
<i>S. pyogenes</i> (group A)	8/7	31/31	38	1	0	2,168	38/38 (100)	2,168/2,169 (99.9)
Total	1,331/1,346	191/189	1,494	28	41	16,093	1,494/1,535 (97.3)	16,093/16,121 (99.8)
Gram-negative bacteria								
<i>A. baumannii</i>	16/14	40/37	51	5	0	2,151	51/51 (100)	2,151/2,156 (99.8)
<i>Enterobacteriaceae</i>	307/310	187/188	490	4	8	1,705	490/498 (98.4)	1,705/1,709 (99.8)
<i>E. cloacae</i> complex	24/22	17/17	38	3	1	2,165	38/39 (97.4)	2,165/2,168 (99.9)
<i>E. coli</i>	149/148	6/5	150	5	3	2,049	150/153 (98.0)	2,049/2,054 (99.8)
<i>K. oxytoca</i>	6/6	54/58	59	1	5	2,142	59/64 (92.2) <sup>c</sup>	2,142/2,143 (99.9)
<i>K. pneumoniae</i>	74/71	37/34	102	9	3	2,093	102/105 (97.1)	2,093/2,102 (99.6)
<i>Proteus</i>	22/22	17/17	39	0	0	2,168	39/39 (100)	2,168/2,168 (100)
<i>S. marcescens</i>	22/22	55/55	76	1	1	2,129	76/77 (98.7)	2,129/2,130 (99.9)
<i>H. influenzae</i>	8/8	35/35	43	0	0	2,164	43/43 (100)	2,164/2,164 (100)
<i>N. meningitidis</i>	1/1	35/35	36	0	0	2,171	36/36 (100)	2,171/2,171 (100)
<i>P. aeruginosa</i>	52/52	0/0	51	1	1	2,154	51/52 (98.1)	2,154/2,155 (99.9)
Total	681/676	483/481	1135	29	22	23,091	1,135/1,157 (98.1)	23,091/23,120 (99.9)

Drug-Resistant Organisms	Isolates detected: BCID/comparator		No. of results: BCID/comparator				Sensitivity or PPA <sup>a</sup> : TP/(TP + FN) (%)	Specificity or NPA <sup>a</sup> : TN/(TN + FP) (%)
	Clinical arm	Seeded arm	TP +/+	FP +/-	FN -/+	TN -/-		
<i>mecA</i> in association with:								
All <i>Staphylococcus</i> isolates detected <sup>b</sup>	491/494	2/2	488	5	8	281	488/496 (98.4)	281/286 (98.3)
<i>Staphylococcus</i> and <i>S. aureus</i> isolates detected	137/139	0/0	137	0	2	118	137/139 (98.6)	118/118 (100)
<i>vanA/B</i> in association with <i>Enterococcus</i> isolates detected	36/36	28/28	64	0	0	67	64/64 (100) <sup>c</sup>	67/67 (100)
<i>bla</i> <sub>KPC</sub> in association with <i>Enterobacteriaceae</i> and/or <i>A. baumannii</i> and/or <i>P. aeruginosa</i> isolates detected	6/6 <sup>d</sup>	33/33	39	0	0	558	39/39 (100) <sup>e</sup>	558/558 (100)

Only KPC!



## Emergence of *Klebsiella pneumoniae* co-producing NDM-1, OXA-48, CTX-M-15, CMY-16, QnrA and ArmA in Switzerland IJAA, 2014

Salome N. Seiffert<sup>a,b,c</sup>, Jonas Marshall<sup>d</sup>, Vincent Perreten<sup>b</sup>, Alessandra Carattoli<sup>a,b,e</sup>, Hansjakob Furrer<sup>d</sup>, Andrea Endimiani<sup>a,\*</sup>



## Real-time long read whole genome sequencing (Point of care, POC)

Antibiotic MIC, in µg/mL (interpretation)<sup>a</sup>


	≥128 (R)
Piperacillin/tazobactam	
Ticarcillin/clavulanic acid	≥256 (R)
Cefoxitin	≥128 (NA)
Cefpodoxime	≥64 (R)
Ceftriaxone	≥256 (R)
Cefotaxime	≥128 (R)
Cefotaxime/clavulanic acid	≥128 (NA)
Ceftazidime	≥256 (R)
Ceftazidime/clavulanic acid	≥256 (NA)
Cefepime	≥32 (R)
Aztreonam	≥32 (R)
Imipenem	4 (R)
Meropenem	≥16 (R)
Ertapenem	≥8 (R)
Doripenem	≥4 (R)
Gentamicin	≥32 (R)
Tobramycin	≥16 (R)
Amikacin	≥64 (R)
Ciprofloxacin	≥4 (R)
Levofloxacin	≥16 (R)
Trimethoprim/sulfamethoxazole	≥8 (R)
Colistin	≤0.125 (S)
Polymyxin B	≤0.125 (NA)
Doxycycline	16 (NA)
Minocycline	8 (NA)
Tigecycline	0.5 (S)
Fosfomycin <sup>b</sup>	2 (S)

Resistance genes (genetic background or characteristics)	bla <sub>SHV-1</sub> (NS) bla <sub>TEM-1</sub> (NS) bla <sub>CTX-M-15</sub> (ISEcp1) bla <sub>CMY-16</sub> (ISEcp1) bla <sub>OXA-48</sub> (Tn1999.2) bla <sub>NDM-1</sub> (ISAb125) OmpK35 (stop codon) armA (NS) GyrA (Ser83Tyr; Asp87Asn) ParC (Ser80Ile) qnrA bla <sub>OXA-7-like</sub> , arr-1, dfrA12, dfrA14, sul1, sul2, cmlA1, floR, aadA1, tet(A), Int11 <sup>c</sup>
Plasmid incompatibility (Inc) group	A/C, R, I/M
MLST	ST101


IncL/M  
OXA-48  
**Carbapenems**

IncR  
CTX-M-15  
**Penicillins  
Cephalosporins  
Monobactams  
Bactrim**

IncA/C  
NDM-1  
CMY-16  
**Penicillins  
Cephalosporins  
Monobactams  
Carbapenems  
Bactrim  
Tetracyclines  
Chloramphenicol  
Aminoglycosides  
Quinolones**



150-600  
Euro



2-3 weeks work for the PhD student

1-2 hrs enough reads/coverage  
To reports all resistance genes

# Conclusions

**Carbapenem-resistant strains have high prevalence and still increase**

- Most of them are carbapenemase producers
- KPC, NDM, and OXA-48 are the most frequent enzymes
- The genetic background and special clones support their spread

**Standard phenotypic tests do not detect all carbapenemase producers**

- The clinical breakpoints do not identify 15-20% of strains
- Screening tests may improve detection (if implemented)
  - They are cheap, but time consuming (1-3 days)

**Non-molecular rapid tests**

- They are cheap and relatively fast, but provide only a generic detection

**Molecular rapid tests**

- Most identify only one or few carbapenemase genes
- Most do not identify species and other important genes (ideal system!)
- Only some are validated for clinical samples
- They are expensive

# THANK YOU!

## Institute for Infectious Diseases - University of Bern, Switzerland

- Dr. med. Baharak Babouee-Flury (PostDoc)
- Odette J. Bernasconi (PhD student)
- Fernanda Pimentel (PhD student)
- Regula Tinguely (Lab Tech)
- Thomas Büdel (Lab Tech)
- Maximilian Scheidegger (Master student)



\* H. Sigrist Foundation (Bern)

