

Vienna – April 22<sup>nd</sup>, 2017

(Hall B - 13:30-15:30)

**Session SY021:**


*Reporting Antibigram Data*

# Direct Reporting of Resistance Mechanisms


**Prof. Andrea Endimiani, MD, PhD**

Institute for Infectious Diseases (IFIK)

University of Bern, Switzerland



Adapted from Nordmann P. *et al.*, EID 18(9), Sept 2012



# MALDI-TOF MS

**Multicenter evaluation of the Sepsityper™ extraction kit and MALDI-TOF MS for direct identification of positive blood culture isolates using the BD BACTEC™ FX and VersaTREK® diagnostic blood culture systems**

Schieffer KM *et al.*, J Appl Micro; 116, 2014



Sepsityper extraction kit  
(Bruker Daltonics)

~15 min

(vs. 7-45 min for lab methods)

~80-85% efficient extraction  
(vs. 85-90% for lab methods)

n=411 +BCs

Accuracy of Sepsityper extraction/MALDI-TOF MS (SE-MALDI) identification method

Organism group	Total
Staphylococci	171 (94.2)
Streptococci	15 (93.3)
Enterococci	25 (92.0)
Other Gram-positive bacteria	14 (64.3)
Enterobacteriaceae	88 (98.9)
Other Gram-negative bacilli	4 (100)
Anaerobic bacteria	2 (0.0)
Yeasts	6 (50.0)
Mixed organisms	33 (93.9)
Total	358 (92.7)

Incl. those with score <1.6


ID in ~1h

Reduction of 6-83 hrs

**Rapid detection and identification of strains carrying carbapenemases directly from positive blood cultures using MALDI-TOF MS**

Hoyos-Mallecot *et al.*, J Micro Meth; 105, 2014

19 carba producers (21 controls)  
Sens and Spec: 100% and 90%



Generic activity in ~4.5 hrs


Can be reduced to ~30 min-2 hrs

# Colorimetric tests

## Rapid Detection of ESBL-Producing *Enterobacteriaceae* in Blood Cultures


Dortet L. *et al.*, Emerg Infect Dis; 21:3, 2015

### NDP test



## Rapid Detection of Polymyxin-Resistant *Enterobacteriaceae* from Blood Cultures

Jayol A. *et al.*, J Clin Micro; 54:9, 2016



**Final results in 2-4 hrs**

**Sens and Spec: 98% and 100%**

## Rapid detection of carbapenemase-producing *Enterobacteriaceae* from blood cultures


Dortet L. *et al.*, Clin Microbiol Inf; 20, 2014

**Preparation from +BCs: ~3 hrs**

**Reaction in 5 min-1 h (max 2 hrs)**

**Sens and Spec 98% and 100%**

**The RAPIDEC Carba NP ®Test**



**RAPID POLYMYXIN NP**

**~10 Euro**



## Gram-negatives

Antibiotics	Mechanism	Most frequent protein/genes
Cephalosporins	ESBLs	CTX-Ms
		TEMs and SHVs (not all are ESBLs)
		PERs, VEBs, GESs, SFO, etc.
	Plasmidic AmpCs	CMYs, DHAs
		FOXs, LATs, BILs, MOXs, MIRs, ACTs
Carbapenems	Chromos. AmpCs	Promoter region ( <i>E. coli</i> )
		AmpD ( <i>Enterobacter</i> spp.)
		Class A: - KPCs, GESs
		Class B: - NDMs, VIMs, IMPs
		- SPMs, GIMs, AIMs, SIMs
		Class D: - OXA-48
		- OXA-48-like
		- OXA-23, -24/40, -58
	Porins	OmpK36 ( <i>K. pneumoniae</i> )
Quinolones	QRDR	OmpF/OmpC ( <i>E. coli</i> )
		OmpF/OmpC ( <i>E. cloacae</i> )
		Substitutions in GyrA and ParC
		QnrA/B/S, QepA, aac(6)-Ib-cr
	PMQR determinants	
Aminoglycosides	AMEs	aac(6)-I-like, ant(3)-like, aph(3)-like
	16S rRNA methyl.	ArmA, RmtA/B/C/D/E, NpmA
Polymyxins	Plasmid mediated	Mcr-1/-2-like
	Chromosomal	mgrB ( <i>K. pneumoniae</i> )
Fosfomycin	Plasmid-mediated	fos genes

## Ideal system:

- Species identification
- Resistance mechanisms
- Kit and automated (easy to use)
- High sensitivity/specificity
- Clinical samples (or at least +BCs)
- Rapid
- Cost-effective

## Gram-positives



Antibiotics	Mechanism	Most frequent genes
B-lactams	PBP2	<i>mecA</i> , <i>mecAv</i> , <i>mecB</i> , <i>mecC</i>
	β-lactamases	<i>blaZ</i>
Glycopeptides		<i>vanA</i> , <i>vanB</i> , <i>vanC/D/E/G/Z</i>
Quinolones	QRDR	Substitutions in GyrA and ParC
Aminoglycosides	AMEs	<i>aac(6)-I-like</i> , <i>aph(2)-I-like</i> .....
Macrolides		<i>ermA/B/C</i> , <i>msr</i> , <i>mefA/E</i> , <i>ereA/B</i> ...
Clindamycin		
Lincosamides		<i>lnuA-D</i> , <i>lsaB-E</i>
Streptogramins		<i>vat</i> , <i>vga</i> , <i>vgb</i> genes
Linezolid		<i>cfr</i>
Mupirocin		<i>mupR</i>
Phenicol		<i>cat</i> genes, <i>fexA</i>
Trimethoprim	Folate inhibitors	<i>dfrA</i> , <i>dfrD</i> , <i>dfrG</i> , <i>dfrK</i>
Tetracyclines		<i>tetK</i> , <i>tetL</i> , <i>tetM/O/Q/S</i> .....
Fusidic acid		<i>fusB</i> , <i>fusC</i> , <i>fusD</i>

# Check-Direct Screening, Check-Points

**Rectal  
swab**  
(Copan, ESwab)




- Rapid preparation
- Time to results (<3 hrs)



**BD MAX**



**Reagent strip**



Check-Direct **ESBL**  
Screen for BD MAX™

CTX-M-1 group  
CTX-M-2 group  
CTX-M-9 group  
SHV-ESBL

Check-Direct **CPE**  
Screen for BD MAX™

KPC  
OXA-48-like  
VIM  
NDM

~30 Euro



Performance of the BD MAX™ instrument with Check-Direct CPE real-time PCR for the detection of carbapenemase genes from rectal swabs, in a setting with endemic dissemination of carbapenemase-producing *Enterobacteriaceae* Antonelli *et al.*, DMID, 2016

Florence, Italy  
557 rectal swabs

## Limit of detection (LOD)

Target	LOD BDMix (CFU/ml)	LOD CPMix (CFU/ml)
<i>bla</i> <sub>KPC</sub>	9×10 <sup>5</sup>	9×10 <sup>3</sup>
<i>bla</i> <sub>VIM</sub>	4.5×10 <sup>2</sup>	4.5×10 <sup>2</sup>
<i>bla</i> <sub>OXA-48</sub>	8.5×10 <sup>1</sup>	8.5×10 <sup>1</sup>
<i>bla</i> <sub>NDM</sub>	7.3×10 <sup>4</sup>	7.3×10 <sup>1</sup>

Culture vs. C-D direct CPE (2 ≠ master mix)

## Summary of the results

		Colonization status	
		Positive (n. 29)	Negative (n. 528)
CHSM and BERM	Positive	24 (82.8 %) <sup>a</sup>	0
	Negative	5 (17.2 %)	528 (100 %)
BDMix 18 (3%) unresolved	Positive	27 (93.1 %)	14 (2.7 %)
	Negative	2 (6.9 %)	496 (93.9 %)
	Unresolved	-	18 (3.4 %)
CPMix 6 (1%) unresolved	Positive	29 (100 %)	15 (2.9 %)
	Negative	0	507 (96.0 %)
	Unresolved	-	6 (1.1 %)

ChromID CARBA SMART (w/wo broth enrichment) →  
[5 samples not detected]

TAT from 18-24 hrs (direct culture)  
or 48 hrs (broth enrichment) to 3 hrs

23 KPC  
5 VIM  
1 OXA-48

## Performance of Check-Direct CPE screening for BD MAX

	Positives for <i>bla</i> <sub>KPC</sub>		Positives for <i>bla</i> <sub>VIM</sub>		Positives for <i>bla</i> <sub>OXA-48</sub>		Positives for <i>bla</i> <sub>NDM</sub>	
	BDMix	CPMix	BDMix	CPMix	BDMix	CPMix	BDMix	CPMix
sensitivity	91,3%	100%	100%	100%	100%	100%	-	-
specificity	98,6%	99,6%	99,3%	97,6%	99,3%	99,6%	100%	100%
PLR	65.21	250	142.86	41.67	- <sup>a</sup>	-	-	-
NLR	0.09	0	0	0	-	-	-	-

# Cepheid GeneXpert



- Add aliquot to elution, vortex, transfer to port S
- Insert cartridge to station (overall, **1 min**)
- Run time (**<1 h**)

## Real-time multiplex PCR

- Smart fluidic system
- Filtering and Sonication (DNA)
- Fluorescent-labeled hybr. probes (6 colors)
- Internal control

### Xpert® Carba-R

KPC  
NDM  
OXA-48-like  
VIM  
IMP-1

### Xpert® vanA

### Xpert® MRSA/SA BC

MSSA  
MRSA  
CoNS

*spa, mecA*  
*SCCmec-orfX*

~50 Euro



# Multisite Evaluation of Cepheid Xpert Carba-R Assay for Detection of Carbapenemase-Producing Organisms in Rectal Swabs JCM, 54:7; 2016

M. Tato,<sup>a</sup> P. Ruiz-Garbajosa,<sup>a</sup> M. Traczewski,<sup>b</sup> A. Dodgson,<sup>c</sup> A. McEwan,<sup>c</sup> R. Humphries,<sup>d</sup> J. Hindler,<sup>d</sup> J. Veltman,<sup>e</sup> H. Wang,<sup>f</sup> R. Cantón<sup>a</sup>

4 centers (2 USA, 1 UK, 1 Spain)  
July 2013 - Feb 2014  
633 samples

50% of them spiked at  $1.1 \times 10^2$  to  $1.2 \times 10^3$  CFU/swab (LOD of Xpert)

## Results (obtained in 32-48 min)

Xpert Carba-R assay result	Clinical specimens (n = 383)	Contrived specimens (n = 250)	All specimens (n = 633)
Positive (single and/or combined targets)	42	107	149 23.5%
IMP-1	0	25	25
VIM	2	24	26
NDM	2	23	25
KPC	13	19	32
OXA-48	20	15	35
VIM + OXA-48	4	0	4
NDM + KPC	1	0	1
IMP-1 + NDM	0	1	1
Negative	341	143	484

## Results by individual targets

Xpert Carba-R assay	Reference method (culture plus sequencing)						Total
	IMP-1	VIM	NDM	KPC	OXA-48	Negative	
IMP-1	26	0	0	0	0	0	26
VIM	0	29	0	0	0	1	30
NDM	0	0	26	0	0	1	27
KPC	0	0	0	29	0	4	33
OXA-48	0	0	0	0	38	1	39
Negative	1	2	0	1	2	3,004	3,010
Total	27	31	26	30	40	3,011	3,165

6 out of 154 specimens not detected by Xpert  
6 cases of multiple genes (4 not detected by ref. method)

## Performance vs. reference method

Xpert Carba-R assay	Reference method (culture plus sequencing)		
	No. positive	No. negative	Total No.
Positive	142	7	149
Negative	6	478	484
Total	148	485	633

## Performance for different targets

Target gene	Sensitivity (% [95% CI])	Specificity (% [95% CI])	PPV (%)	NPV (%)
IMP-1	96.3 (81.0–99.9)	100 (99.4–100)	100	99.8
VIM	93.5 (78.6–99.2)	99.8 (99.1–100)	96.7	99.7
NDM	100 (86.8–100)	99.8 (99.1–100)	96.3	100
KPC	96.7 (82.8–99.9)	99.3 (98.3–99.8)	87.9	99.8
OXA-48	95.0 (83.1–99.4)	99.8 (99.1–100)	97.4	99.7

Sensitivity: 96.6%  
Specificity: 98.6%

Comparison of the Next-Generation Xpert MRSA/SA BC Assay and the GeneOhm StaphSR Assay to Routine Culture for Identification of *Staphylococcus aureus* and Methicillin-Resistant *S. aureus* in Positive-Blood-Culture Broths Buchan BW *et al.*, J Clin Microb, 53:3; 2015

Multicenter study (8 hospitals, USA)  
795 BCs (30%, *S. aureus*; 13% MRSA)

Initial test validity

- Xpert: 96.1%
- GeneOhm: 96.3%

Retest:

- Xpert: 99.6%
- GeneOhm: 98.1%

## Identification of *Staphylococcus aureus*

Test	Site	Total no. of specimens tested	No. with result <sup>a</sup> :				Sensitivity (% [95% CI]) <sup>b</sup>	Specificity (% [95% CI])
			TP	FP	TN	FN		
Xpert MRSA/SA BC	A	63	19	0	44	0	100 (82.4–100)	100 (92.0–100)
	B	91	32	0	59	0	100 (89.1–100)	100 (93.9–100)
	C	44	17	0	27	0	100 (80.5–100)	100 (87.2–100)
	D	70	26	0	44	0	100 (86.8–100)	100 (92.0–100)
	E	130	27	2 <sup>c</sup>	100	1 <sup>d</sup>	96.4 (79.7–99.8)	98.0 (92.4–99.6)
	F	211	65	1 <sup>e</sup>	145	0	100 (94.5–100)	99.3 (96.2–100)
	G	126	31	0	95	0	100 (88.8–100)	100 (96.2–100)
	H	57	18	0	39	0	100 (81.5–100)	100 (91.0–100)
	Total	792	235	3	553	1	99.6 (97.7–99.9)	99.5 (98.4–99.9)
GeneOhm StaphSR	Total	782	234	19	527	2	99.2 (97.0–99.9)	96.5 (94.6–97.9)

} Statistically equivalent

## Identification of MRSA


Test	Site	Total no. of specimens tested	No. with result <sup>a</sup> :				Sensitivity (% [95% CI]) <sup>b</sup>	Specificity (% [95% CI])
			TP	FP	TN	FN		
Xpert MRSA/SA BC	A	63	11	0	52	0	100 (71.5–100)	100 (93.2–100)
	B	91	16	0	75	0	100 (79.4–100)	100 (95.2–100)
	C	44	7	0	36	1	87.5 (47.3–99.7)	100 (90.3–100)
	D	70	11	1 <sup>c</sup>	58	0	100 (71.5–100)	98.3 (90.9–100)
	E	130	10	1 <sup>d</sup>	119	0	100 (69.2–100)	99.2 (95.4–100)
	F	211	25	1 <sup>c</sup>	185	0	100 (86.3–100)	99.5 (97.0–100)
	G	126	19	0	106	1	95.0 (73.0–99.7)	100 (96.2–100)
	H	57	4	0	53	0	100 (81.5–100)	100 (91.6–100)
	Total	792	103	3	684	2 <sup>e</sup>	98.1 (93.3–99.8)	99.6 (98.7–99.9)
GeneOhm StaphSR	Total	782	99	15	662	6 <sup>i</sup>	94.3 (88.0–97.9)	97.8 (96.4–98.8)

} Statistically equivalent

# Verigene, Nanosphere



- Load cartridge, consumables, and sample (**5 min**)
- Automated sample preparation and processing
- Place slide from cartridge in reader (**2.5 hrs**)



Microarray approach by using Au-nanoprobe as reporter and silver reduction to enhance signal

## Gram-negatives cartridge

Species	Genus	Resistance
<i>Escherichia coli</i> *	<i>Acinetobacter</i> spp.	CTX-M (ESBL)
<i>Klebsiella pneumoniae</i>	<i>Citrobacter</i> spp.	IMP (carbapenemase)
<i>Klebsiella oxytoca</i>	<i>Enterobacter</i> spp.	KPC (carbapenemase)
<i>Pseudomonas aeruginosa</i>	<i>Proteus</i> spp.	NDM (carbapenemase)
<i>Serratia marcescens</i>		OXA (carbapenemase)
		VIM (carbapenemase)

## Gram-positives cartridge

Species	Genus	Resistance
<i>Staphylococcus aureus</i>	<i>Staphylococcus</i> spp.	<i>mecA</i> (methicillin)
<i>Staphylococcus epidermidis</i>	<i>Streptococcus</i> spp.	<i>vanA</i> (vancomycin)
<i>Staphylococcus lugdunensis</i>	<i>Micrococcus</i> spp.	<i>vanB</i> (vancomycin)
<i>Streptococcus anginosus</i> Group	<i>Listeria</i> spp.	
<i>Streptococcus agalactiae</i>		
<i>Streptococcus pneumoniae</i>		
<i>Streptococcus pyogenes</i>		
<i>Enterococcus faecalis</i>		
<i>Enterococcus faecium</i>		

~50 Euro

# Performance Evaluation of the Verigene Gram-Positive and Gram-Negative Blood Culture Test for Direct Identification of Bacteria and Their Resistance Determinants from Positive Blood Cultures in Hong Kong

Gilman K.H. Siu *et al.*, Plos One, Oct 2015

Multicenter study (4 hospitals)

364 BCs (114 Gram-pos; 250 Gram-neg)

Jan 2014 - May 2014

Initial test validity: 95.6%

## ID for Gram-positives (agreement 89.6%)

Organisms	No. (%) of isolates					Sensitivity (%)	Specificity (%)
	Total	Correctly identified	Not detected	Misidentified	No Call		
<i>Staphylococcus</i> spp.	73 (60.3)	69 (94.5)	-	-	4	100	100
<i>S. aureus</i>	48 (39.7)	47 (97.9)	-	-	1	100	100
<i>S. epidermidis</i>	4 (3.3)	4 (100)	-	-	-	100	100
<i>Streptococcus</i> spp.	27 (22.3)	25 (92.6)	-	-	2	100	98.9
<i>S. pyogenes</i>	2 (1.7)	2 (100)	-	-	-	100	100
<i>S. agalactiae</i>	3 (2.5)	3 (100)	-	-	-	100	100
<i>S. anginosus</i> gr.	1 (0.8)	1 (100)	-	-	-	100	100
<i>S. pneumoniae</i>	7 (5.8)	7 (100)	-	-	-	100	99.1
<i>E. faecalis</i>	7 (5.8)	4 (57.1)	2 (28.6)	-	1	66.7	100
<i>E. faecium</i>	5 (4.1)	2 (40)	2 (40)	-	1	50	100
<i>Listeria</i>	1 (0.8)	1 (100)	-	-	-	100	100

## ID for Gram-negatives (agreement 90.5%)

Organisms	No. (%) of isolates					Sensitivity (%)	Specificity (%)
	Total	Correctly identified	Not detected	Misidentified	No Call		
<i>E. coli</i>	165 (59.1)	158 (95.8)	6 (3.6)	1 <sup>c</sup>	-	95.8	100
<i>K. pneumoniae</i>	52 (18.6)	36 (69.2)	14 (27) <sup>d</sup>	2 <sup>e</sup>	-	69.2 <sup>d</sup>	100
<i>P. aeruginosa</i>	17 (6.1)	13 (76.5)	3 (17.6)	-	1	81.3	100
<i>Proteus</i> spp.	10 (3.6)	9 (90)	-	-	1	100	100
<i>Enterobacter</i> spp.	10 (3.6)	8 (80)	-	-	2	100	99.3
<i>Acinetobacter</i> spp.	4 (1.4)	4 (100)	-	-	-	100	100
<i>K. oxytoca</i>	3 (1.1)	-	3 (100)	-	-	0	99.3
<i>Citrobacter</i> spp.	1 (0.4)	1 (100)	-	-	-	100	100
<i>Serratia marcescens</i>	1 (0.4)	1 (100)	-	-	-	100	100

*A. hydrophila*, *M. morganii*, *Salmonella* spp., *A. faecalis*, *B. pseudomallei*, *H. influenzae*, *R. planticola*, *P. putida*, *S. maltophilia* not detectable with the system (n=16; 5.5%)

## Drug-Resistant Organisms

Drug resistant Organisms	No. of isolates					Sensitivity	Specificity (%)
	Total	Correctly Detected	Not Detected	No Call			
<b>Gram Positive</b>							
MRSA	27	26	0	1	100	100	
MRSE	4	4	0	0	100	100	
VRE	1	1	0	0	100	100	
<b>Gram Negative</b>							
Cefotaxime resistant <i>Enterobacteriaceae</i> (including ESBL producers)	61	38	22	1	63.3	100	
ESBL producing <i>Enterobacteriaceae</i>	46	38	7	1	84.4	100	
MDR <i>Acinetobacter</i>	3	3	0	0	100	100	
Carbapenem-resistant <i>Pseudomonas</i>	2	0	2	0	0	100	
<b>Total</b>	<b>98</b>	<b>72</b>	<b>24</b>	<b>2</b>	<b>75.0</b>	<b>100</b>	

## Time to Results (40-99 hrs faster than routine)

Organisms	No. of Isolates <sup>a</sup>	Δ Time to Result <sup>b</sup>				p-value
		Average Time to Identification by Culture-Based Method (h)	Average Time to Result by Verigene Test (h)	Average (h)	Range (h)	
<i>Staphylococcus</i> spp.	23	63.76	2.35	61.41	22.65-494.65	p<0.001
MSSA	7	49	2.35	46.65	22.65-79.65	p = 0.012
MRSA	10	93.45	2.35	91.1	23.65-494.65	p<0.001
CNS	8	57.19	2.35	54.84	29.65-97.65	p<0.001
<i>Streptococcus</i> spp.	13	63.86	2.35	61.51	22.65-128.15	p<0.001
<i>S. pneumoniae</i>	3	42.83	2.35	40.48	22.65-56.65	p = 0.017
β-haemolytic Strept.	4	71.62	2.35	69.27	33.65-128.15	p<0.001
Viridians Group	6	55.08	2.35	52.73	25.65-83.65	p<0.001
<i>Enterococcus</i> spp.	2	78.5	2.35	76.15	53.15-99.15	p<0.001
VSE	1	55.5	2.35	53.15	n/a	n/a
VRE	1	101.5	2.35	99.15	n/a	n/a
<i>Enterobacteriaceae</i>	79	45.98	1.88	44.1	18.12-104.12	p<0.001
<i>P. aeruginosa</i>	7	52.51	1.88	50.63	28.62-94.62	p<0.001
<i>Acinetobacter</i> spp.	1	76	1.88	74.12	n/a	n/a

# Clinical Impact of Laboratory Implementation of Verigene BC-GN Microarray-Based Assay for Detection of Gram-Negative Bacteria in Positive Blood Cultures

Walker T. *et al.*, J Clin Microbiol; 54:7, 2016

Los Angeles  
May 2013 - Nov 2013  
Dec 2014 - May 2015

6-month period before Verigene (n=98) vs. 6-month period after Verigene (n=97)

## Clinical outcome

	Pre-BC-GN	Post-BC-GN	P value
Mean time from initial Gram stain to BC-GN identification, h	NA <sup>a</sup>	3.5	NA
Mean time from initial Gram stain to organism identification, h	37.9	10.9	<0.001 <sup>b</sup>
Mean time from initial Gram stain to effective therapy, h			
All cases	10.2	6.5	0.12 <sup>b</sup>
Cases on suboptimal empirical therapy	30.3	19.1	0.12 <sup>b</sup>
No. of cases in which therapy was de-escalated	33	37	0.66 <sup>c</sup>
Mean time from initial Gram stain to de-escalation, h	40.9	34.1	0.14 <sup>b</sup>
Recurrence of bacteremia, no. (%)	8 (8.2)	3 (3.1)	0.21 <sup>c</sup>
Mean total length of stay in hospital, days	15.2	18.0	0.52 <sup>b</sup>
Mean length of hospital stay after positive culture, days	9.7	9.4	0.87 <sup>b</sup>
Mean length of stay in ICU, days	16.2	12.0	0.03 <sup>b</sup>
30-day mortality, no. (%)	19 (19.2)	8 (8.1)	0.04 <sup>c</sup>
ESBL cases, no.	15	11	0.53 <sup>c</sup>
Length of stay in hospital, days	12.0	13.5	0.59 <sup>b</sup>
Mean time to effective therapy, h (no.) <sup>d</sup>	41.4 (9)	7.3 (9)	0.04 <sup>b</sup>
30-day mortality, no. (%)	4 (26.7)	0 (0)	0.11 <sup>c</sup>

Saved 11,661 USD/case ←

## Factors associated to length of ICU stay

Factor	Univariate analysis <sup>b</sup>			Multivariate analysis <sup>b</sup>		
	HR	95% CI	P value	HR	95% CI	P value
Intervention	1.62	1.03–2.54	0.04	1.79	1.07–2.98	0.03
Preinfection length of stay	0.97	0.95–0.99	0.001	0.97	0.95–0.98	0.0001
Age	1.01	0.99–1.03	0.22	1.0	0.98–1.01	0.90
Sex (female)	0.64	0.40–1.03	0.06	0.69	0.41–1.15	0.15
Time to effective therapy	1.01	0.78–1.32	0.94	1.02	0.75–1.38	0.90
Disease severity (APACHE II)	1.00	0.97–1.02	0.74	0.97	0.93–1.0	0.06
Mortality	1.03	0.57–1.65	0.91	1.43	0.83–2.48	0.20

## Multiv. log. regression for association with 30-d mortality




Factor	Odds ratio	95% CI	P value
Intervention	0.81	0.67–0.98	0.03
Length of stay, ICU	1.00	1.00–1.00	0.94
Age	1.01	1.00–1.01	0.11
Sex (female)	1.07	0.88–1.30	0.48
APACHE II	1.01	0.99–1.02	0.36



# BioFire FilmArray



- Preparation of the pouch
- Add pouch to FilmArray station (overall, **2 min**)
- Run time of about **1 h**



FilmArray® BCID Panel targets:

~100 Euro

Gram+ Bacteria	Gram- Bacteria
<i>Enterococcus</i> <i>Listeria monocytogenes</i> <i>Staphylococcus</i> <i>Staphylococcus aureus</i> <i>Streptococcus</i> <i>Streptococcus agalactiae</i> <i>Streptococcus pyogenes</i> <i>Streptococcus pneumoniae</i>	<i>Acinetobacter baumannii</i> <i>Haemophilus influenzae</i> <i>Neisseria meningitidis</i> <i>Pseudomonas aeruginosa</i> <i>Enterobacteriaceae</i> <i>Enterobacter cloacae</i> complex <i>Escherichia coli</i> <i>Klebsiella oxytoca</i> <i>Klebsiella pneumoniae</i> <i>Proteus</i> <i>Serratia marcescens</i>
Yeast	Antibiotic Resistance
<i>Candida albicans</i> <i>Candida glabrata</i> <i>Candida krusei</i> <i>Candida parapsilosis</i> <i>Candida tropicalis</i>	mecA - methicillin resistance vanA/B - vancomycin resistance KPC - carbapenem resistance



## Species ID

Species ID	Isolates detected <sup>a</sup> :		No. of results:				Sensitivity or PPA <sup>b</sup> : TP/(TP + FN) (%)	Specificity or NPA <sup>b</sup> : TN/(TN + FP) (%)
	BCID/comparator		BCID/comparator					
	Clinical arm	Seeded arm	TP +/+	FP +/-	FN -/+	TN -/-		
Gram-positive bacteria								
<i>Enterococcus</i>	102/101	29/29	127	4	3	2,073	127/130 (97.7)	2,073/2,077 (99.8)
<i>L. monocytogenes</i>	0/0	36/36	36	0	0	2,171	36/36 (100)	2,171/2,171 (100)
<i>Staphylococcus</i>	780/797	2/1	770	12	28	1,397	770/798 (96.5)	1,397/1,409 (99.1)
<i>S. aureus</i>	257/257	0/0	253	4	4	1,946	253/257 (98.4)	1,946/1,950 (99.8)
<i>Streptococcus</i>	140/141	63/62	198	5	5	1,999	198/203 (97.5)	1,999/2,004 (99.8)
<i>S. agalactiae</i> (group B)	18/18	18/18	36	0	0	2,171	36/36 (100)	2,171/2,171 (100)
<i>S. pneumoniae</i>	26/25	12/12	36	2	1	2,168	36/37 (97.3)	2,168/2,170 (99.9)
<i>S. pyogenes</i> (group A)	8/7	31/31	38	1	0	2,168	38/38 (100)	2,168/2,169 (99.9)
Total	1,331/1,346	191/189	1,494	28	41	16,093	1,494/1,535 (97.3)	16,093/16,121 (99.8)
Gram-negative bacteria								
<i>A. baumannii</i>	16/14	40/37	51	5	0	2,151	51/51 (100)	2,151/2,156 (99.8)
<i>Enterobacteriaceae</i>	307/310	187/188	490	4	8	1,705	490/498 (98.4)	1,705/1,709 (99.8)
<i>E. cloacae</i> complex	24/22	17/17	38	3	1	2,165	38/39 (97.4)	2,165/2,168 (99.9)
<i>E. coli</i>	149/148	6/5	150	5	3	2,049	150/153 (98.0)	2,049/2,054 (99.8)
<i>K. oxytoca</i>	6/6	54/58	59	1	5	2,142	59/64 (92.2) <sup>c</sup> 98%	2,142/2,143 (99.9)
<i>K. pneumoniae</i>	74/71	37/34	102	9	3	2,093	102/105 (97.1)	2,093/2,102 (99.6)
<i>Proteus</i>	22/22	17/17	39	0	0	2,168	39/39 (100)	2,168/2,168 (100)
<i>S. marcescens</i>	22/22	55/55	76	1	1	2,129	76/77 (98.7)	2,129/2,130 (99.9)
<i>H. influenzae</i>	8/8	35/35	43	0	0	2,164	43/43 (100)	2,164/2,164 (100)
<i>N. meningitidis</i>	1/1	35/35	36	0	0	2,171	36/36 (100)	2,171/2,171 (100)
<i>P. aeruginosa</i>	52/52	0/0	51	1	1	2,154	51/52 (98.1)	2,154/2,155 (99.9)
Total	681/676	483/481	1135	29	22	23,091	1,135/1,157 (98.1)	23,091/23,120 (99.9)

## Evaluation of the FilmArray Blood Culture Identification Panel: Results of a Multicenter Controlled Trial

Salimnia H.. *et al.*, J Clin Microbiol; 54:3, 2016

8 centers in USA  
July 2012 - Feb 2014  
2,207 BC samples

- 1,568 clinical  
- 639 seeded

FilmArray performed:

- BC+, within 8 hrs

- 52% clinical samples  
- 66% seeded cultures

- Aliquot frozen (<8 hrs)

Only 1.9% needed repeat testing

## Drug-Resistant Organisms

gene(s)	clinical arm	seeded arm	TP +/+	FP +/-	FN -/+	TN -/-	TP/(TP + FN) (%)	TN/(TN + FP) (%)
<b><i>mecA</i> in association with:</b>								
All <i>Staphylococcus</i> isolates detected <sup>b</sup>	491/494	2/2	488	5	8	281	488/496 (98.4)	281/286 (98.3)
<i>Staphylococcus</i> and <i>S. aureus</i> isolates detected	137/139	0/0	137	0	2	118	137/139 (98.6)	118/118 (100)
<b><i>vanA/B</i> in association with</b>								
<i>Enterococcus</i> isolates detected	36/36	28/28	64	0	0	67	64/64 (100) <sup>c</sup>	67/67 (100)
<b><i>bla</i><sub>KPC</sub> in association with</b>								
<i>Enterobacteriaceae</i> and/ or <i>A. baumannii</i> and/ or <i>P. aeruginosa</i> isolates detected	6/6 <sup>d</sup>	33/33	39	0	0	558	39/39 (100) <sup>e</sup>	558/558 (100)

# Clinical and economic impact of antimicrobial stewardship interventions with the FilmArray blood culture identification panel *Pardo J et al., DMID; 84, 2016*


Gainesville, Florida

Intervention group (FilmArray)  
(Aug 2013 - Jan 2014; n=84)

vs.

Historical control group  
(Jan 2012 - Jun 2013; n=252)

55% CoNS



Time to discharge after  
CoNS-contaminated BCs

Variable	All Patients (N = 336)			Discharged Within 6 Days of CoNS Contamination (N = 100)		
	Control group (N = 252)	BCID group (N = 84)	P	Control group (N = 75)	BCID group (N = 25)	P
Median LOS (days [IQR])	7.9 (3.8–16.3)	7.4 (3–11.7)	0.178	2.9 (2.1–4.4)	2.3 (1.5–3.1)	0.008
Median ICU Days (IQR)	0 (0–5)	0 (0–1)	0.026	0 (0–0)	0 (0–0)	0.838
In-hospital mortality (No. [%])	37 (15)	5 (6)	0.036	4 (6)	0 (0)	0.239
Median Costs (\$ [IQR])						
Total	15324 (5517–37305)	12241 (3080–21027)	0.030	3370 (1537–5190)	1645 (585–2952)	0.016
Pharmacy	2119 (598–5503)	1538 (366–3824)	0.073	429 (157–712)	214 (75–361)	0.010
Ward room	2806 (0–7427)	3972 (947–8041)	0.514	1852 (891–2778)	947 (0–1828)	0.025
ICU room	1222 (0–10218)	0 (0–3756)	0.008	0 (0–0)	0 (0–0)	0.659
Laboratory	708 (241–1804)	656 (254–1445)	0.692	160 (83–253)	204 (160–321)	0.020

129 USD

## Unnecessary vancomycin use

Variable	MSSA (N = 32)			CoNS Contaminants (N = 184)		
	Control Group (N = 24)	BCID Group (N = 8)	P	Control Group (N = 138)	BCID Group (N = 46)	P
Vancomycin started (No. [%])	24 (100)	8 (100)	>0.99	100 (72)	32 (70)	0.705
1 dose only (No. [%])	0/24 (0)	4/8 (50)	0.002	12/100 (12)	7/32 (22)	0.133
Median duration (h [IQR])	66.0 (54.1–78.1)	14.0 (1–20)	<0.001	66.9 (37.6–135)	39.3 (14.7–96.5)	0.100

# Amplex Eazyplex




- Preparation (**5 min**)
- No DNA extraction
- Run time (**<30 min**)

## LAMP

Loop-mediated isothermal Amplification

Real-time fluorescent measurement



~55 Euro

### eazyplex® MRSA

*S. aureus*  
mecA  
mecC  
*S. epidermidis*

### eazyplex® VRE

### eazyplex® SuperBug mcr-1

### eazyplex® SuperBug CRE

KPC,  
NDM,  
OXA-48  
and OXA 181,  
VIM, as well as  
CTX-M-1 and  
CTX-M-9 group

### eazyplex® SuperBug complete

#### Version A

KPC  
NDM  
OXA-48  
VIM  
OXA-23 group  
OXA-40 group  
**OXA-58 group**

#### Version B

KPC  
NDM  
OXA-48  
VIM  
OXA-23 group  
OXA-40 group  
**OXA-181**

# Extended-spectrum $\beta$ -lactamase (ESBL) detection directly from urine samples with the rapid isothermal amplification-based eazyplex® SuperBug CRE assay: Proof of concept☆ Hinic V. *et al.*, J Microb Meth; 119, 2015

50 urine samples

33 ESBL+ and/or carbapenemases

- Vitek 2 and Etest ESBL/AmpC strips
- Eazyplex on colonies

No	Species (phenotype) isolated from urine	Quantity in urine culture (CFU/ml)	eazyplex® SuperBug CRE			
			Bacterial isolate		Urine	
			Target detected	Threshold time (min:sec)	Target detected	Threshold time (min:sec)
1	<i>Escherichia coli</i> (ESBL)	$10^6$	CTX-M-1 group	04:30	CTX-M-1 group	05:45
2	<i>Escherichia coli</i> (ESBL)	$10^6$	CTX-M-1 group	04:45	CTX-M-1 group	05:30
3	<i>Escherichia coli</i> (ESBL)	$10^5$	CTX-M-1 group	04:30	CTX-M-1 group	08:45
4	<i>Escherichia coli</i> (ESBL)	$10^4$	CTX-M-1 group	04:45	CTX-M-1 group	09:00
5	<i>Escherichia coli</i> (ESBL), <i>Enterococcus faecalis</i>	$10^6, 10^5$	CTX-M-1 group	05:00	CTX-M-1 group	05:15
6	<i>Escherichia coli</i> (ESBL)	$10^4$	CTX-M-1 group	04:30	CTX-M-1 group	10:45
7	<i>Escherichia coli</i> (ESBL), <i>Escherichia coli</i> (no $\beta$ -lactamase)	$10^4, 10^4$	CTX-M-1 group	04:54	CTX-M-1 group	08:45
8	<i>Escherichia coli</i> (ESBL)	$10^6$	CTX-M-1 group	05:15	CTX-M-1 group	07:30
9	<i>Escherichia coli</i> (ESBL)	$10^6$	CTX-M-1 group	04:45	CTX-M-1 group	07:00
10	<i>Escherichia coli</i> (ESBL), <i>Escherichia coli</i> (no $\beta$ -lactamase)	$10^6, 10^5$	CTX-M-1 group	04:30	CTX-M-1 group	06:45
11	<i>Escherichia coli</i> (ESBL), <i>Enterococcus faecium</i>	$10^4, 10^5$	CTX-M-1 group	04:45	CTX-M-1 group	11:00
12	<i>Escherichia coli</i> (ESBL)	$10^5$	CTX-M-1 group	04:30	CTX-M-1 group	07:30
13	<i>Escherichia coli</i> (ESBL), gram-positive flora	$10^5, 10^4$	CTX-M-1 group	04:30	CTX-M-1 group	07:45
14	<i>Escherichia coli</i> (ESBL)	$10^6$	CTX-M-1 group	04:15	CTX-M-1 group	05:30
15	<i>Escherichia coli</i> (ESBL)	$10^5$	CTX-M-1 group	05:15	CTX-M-1 group	06:00
16	<i>Escherichia coli</i> (ESBL)	$10^6$	CTX-M-1 group	05:00	CTX-M-1 group	05:15
17	<i>Escherichia coli</i> (ESBL), <i>Klebsiella pneumoniae</i> (no $\beta$ -lactamase)	$10^4, 10^3$	CTX-M-1 group	05:15	CTX-M-1 group	06:00
18	<i>Escherichia coli</i> (ESBL), <i>Enterococcus faecalis</i>	$10^4, 10^4$	CTX-M-1 group	05:30	CTX-M-1 group	14:00
19	<i>Escherichia coli</i> (ESBL)	$10^6$	CTX-M-1 group	04:15	CTX-M-1 group	04:30
20	<i>Klebsiella pneumoniae</i> (ESBL)	$10^6$	CTX-M-1 group	04:30	CTX-M-1 group	04:30
21	<i>Klebsiella pneumoniae</i> (ESBL)	$10^6$	CTX-M-1 group	05:00	CTX-M-1 group	06:30
22	<i>Klebsiella pneumoniae</i> (ESBL)	$10^4$	CTX-M-1 group	04:45	CTX-M-1 group	08:30
23	<i>Escherichia coli</i> (ESBL)	$10^6$	CTX-M-9 group	05:45	CTX-M-9 group	06:30
24	<i>Escherichia coli</i> (ESBL)	$10^5$	CTX-M-9 group	06:00	CTX-M-9 group	12:45
25	<i>Escherichia coli</i> (ESBL), gram-positive flora	$10^6, 10^5$	CTX-M-9 group	05:15	CTX-M-9 group	05:45
26	<i>Escherichia coli</i> (ESBL)	$10^5$	CTX-M-9 group	04:45	CTX-M-9 group	08:15
27	<i>Escherichia coli</i> (ESBL)	$10^4$	CTX-M-9 group	05:30	CTX-M-9 group	12:30
28	<i>Escherichia coli</i> (ESBL)	$10^6$	CTX-M-9 group	07:30	CTX-M-9 group	11:15
29	<i>Escherichia coli</i> (ESBL), mixed flora 2 x	$10^5, 10^3$	CTX-M-9 group	08:15	CTX-M-9 group	13:45
30	<i>Klebsiella pneumoniae</i> (ESBL)	$10^6$	CTX-M-9 group	06:15	CTX-M-9 group	07:30
31	<i>Escherichia coli</i> (AmpC)	$10^5$	negative		negative	
32	<i>Escherichia coli</i> (AmpC), gram-positive flora 2 x	$10^5, 10^4$	negative		negative	
33	<i>Escherichia coli</i> (AmpC)	$10^6$	negative		negative	
34	<i>Escherichia coli</i> (AmpC), gram-positive flora	$10^4, 10^5$	negative		negative	
35	<i>Enterobacter cloacae</i> (AmpC)	$10^5$	negative		negative	
36	Gram-positive flora 4 x	$10^5$	negative		negative	
37	<i>Enterococcus faecalis</i>	$10^5$	negative		negative	
38	<i>Citrobacter koseri</i> (no $\beta$ -lactamase), <i>Aerococcus urinae</i>	$10^6, 10^4$	negative		negative	
39	<i>Lactobacillus crispatus</i> , gram-positive flora 3 x	$10^5, 10^4$	negative		negative	
40	<i>Candida glabrata</i>	$10^5$	negative		negative	
41	<i>Escherichia coli</i> (no $\beta$ -lactamase)	$10^5$			KPC, NDM, VIM, OXA-48, OXA-181, CTX-M-1 group, CTX-M-9 group <sup>a</sup>	16:30, 17:00, 16:45, 18:30, 19:15, 16:00, 19:15
42	<i>Pseudomonas aeruginosa</i> , <i>Morganella morganii</i> (AmpC)	$10^5, 10^6$			negative	
43	<i>Klebsiella pneumoniae</i> (no $\beta$ -lactamase)	$10^6$			negative	
44	<i>Proteus mirabilis</i> (no $\beta$ -lactamase)	$10^5$			negative	
45	<i>Enterococcus faecalis</i> , <i>Citrobacter freundii</i> -group (AmpC), <i>Pseudomonas aeruginosa</i> , <i>Achromobacter xylosoxidans</i>	$10^5, 10^5, 10^4, 10^5$			negative	
46	<i>Escherichia coli</i> (no $\beta$ -lactamase)	$10^6$			negative	
47	<i>Escherichia coli</i> (no $\beta$ -lactamase)	$10^6$			KPC, NDM, VIM, OXA-48, OXA-181,	08:15, 08:30, 08:30, 08:30, 09:30, 08:00,

## LOD in Urine

Determination of analytical sensitivity of direct urine testing by using a CTX-M-1-group-producing *E. coli* strain (701541/14).

No. of CFU/ml	Threshold time (min:sec)	Result
$10^5$	08:45	positive
$10^4$	11:00	positive
$10^3$	16:00	positive
$10^2$	0	negative

First run with 50 urines:

- 2 invalid results
- 1 false-positive NDM (19 min)

Sensitivity: 100%  
Specificity: 97.9%

Avg. 8 min

# Evaluation of loop-mediated isothermal amplification for the rapid identification of bacteria and resistance determinants in positive blood cultures

Rödel J. *et al.*, Eur J Clin Microbiol Infect Dis; Jan 2017

Jena, Germany

April to Aug 2015

370 BCs + (of which 140 G+)

eazyplex® MRSA

Table 1 Analytes detected by the eazyplex® assays used in this study<sup>a</sup>

MRSA	GP	GN
<i>S. aureus</i>	<i>E. faecalis</i>	<i>E. coli</i>
<i>S. epidermidis</i>	<i>Enterococcus</i> spp.	<i>K. pneumoniae</i>
<i>mecA</i>	<i>S. pneumoniae</i>	Enterobacteriaceae
<i>mecC</i>	<i>Streptococcus</i> spp.	<i>Pseudomonas</i> spp.
	<i>vanA</i>	<i>bla</i> <sub>CTX-M-1</sub> group
	<i>vanB</i>	<i>bla</i> <sub>CTX-M-9</sub> group

In-house LAMP assays (not available)

ID & Resistance	True-positive (n)	True-negative (n)	False-positive (n)	False-negative (n)	No. positive/no. tested (%) (95% CI <sup>a</sup> )	
					Sensitivity	Specificity
<b>MRSA</b>						
<i>S. aureus</i>	31	106	2	0	100 (88.8–100)	98.2 (93.5–99.8)
MRSA <sup>d</sup>	6	25	0	0	100 (54.1–100)	100 (86.3–100)
<i>S. epidermidis</i>	72	60	1	6	92.3 (84–97.1)	98.4 (91.2–100)
<i>mecA</i> (CoNS)	73	30	0	3	96 (88.9–99.2)	100 (88.4–100)
<b>GP</b>						
<i>E. faecalis</i>	12	58	1	0	100 (73.5–100)	98.3 (90.9–100)
<i>Enterococcus</i> spp.	41	27	2	1	97.6 (87.4–99.9)	93.1 (77.2–99.1)
VRE <sup>e</sup>	3	39	0	0	100 (29.2–100)	100 (91–100)
<i>S. pneumoniae</i>	2	69	0	0	100 (15.8–100)	100 (94.8–100)
<i>Streptococcus</i> spp.	20	50	0	1	95.2 (76.2–99.9)	100 (92.9–100)
<b>GN</b>						
<i>E. coli</i>	97	58	2	0	100 (96.3–100)	96.7 (88.5–99.6)
<i>E. coli</i> CTX-M <sup>f</sup>	17	80	0	0	100 (89.5–100)	100 (95.5–100)
<i>E. coli</i> ESBL <sup>g</sup>	17	78	0	2	89.5 (66.9–98.7)	100 (95.3–100)
<i>K. pneumoniae</i>	11	143	2	1	91.7 (61.5–99.8)	98.6 (95.1–99.8)
<i>K. pneumoniae</i> ESBL <sup>h</sup>	4	7	0	0	100 (39.8–100)	100 (59–100)
Enterobacteriaceae	134	12	11	2	98.5 (94.8–99.8)	52.2 (30.6–73.2)
<i>Pseudomonas</i> spp.	9	148	0	0	100 (66.4–100)	100 (97.5–100)

Time to results	Threshold time [min; mean values (SD)]	
	True-positive	False-positive
<b>MRSA</b>		
<i>S. aureus</i>	7 (1.5)	17
<i>mecA</i> ( <i>S. aureus</i> )	11.5 (4)	–
<i>S. epidermidis</i>	11.75 (2.75)	18.75
<i>mecA</i> (CoNS)	9.75 (2.5)	–
IC	7.5 (1)	–
<b>GP</b>		
<i>E. faecalis</i>	6.75 (1)	17.5
<i>Enterococcus</i> spp.	10 (4)	17.5
<i>vanA</i>	7 (1)	–
<i>vanB</i>	–	–
<i>S. pneumoniae</i>	7	–
<i>Streptococcus</i> spp.	13.5 (3.25)	–
IC	9.5 (0.75)	–
<b>GN</b>		
<i>E. coli</i>	7.25 (1.75)	15
<i>K. pneumoniae</i>	13.75 (2.5)	10.75
<i>bla</i> <sub>CTX-M-1</sub> group <sup>a</sup>	8.5 (4)	–
<i>bla</i> <sub>CTX-M-9</sub> group <sup>b</sup>	6.25 (1)	–
Enterobacteriaceae	6.25 (1.75)	15.75 (2.5)
<i>Pseudomonas</i> spp.	8.25 (1.5)	–
IC	9 (1.75)	–


In principle, from +BC ~30 min



# Curetis Unyvero



- Sample transferred to Unyvero tube (2 min)
- Lysis (30 min), transfer to cartridge, plus master mix (3 min)
- Overall, time to results (~5 hrs)



DNA purification, PCR set up, multiplex PCRs with array detection

## Pneumonia panel (P55)

GROUP	PATHOGEN	GENE	RESISTANCE AGAINST
Gram-positive bacteria	<i>Staphylococcus aureus</i>	<i>ctx-M</i>	3rd generation Cephalosporins
	<i>Streptococcus pneumoniae</i>	<i>ermB</i>	Macrolide/Lincosamide
Enterobacteriaceae	<i>Citrobacter freundii</i>	<i>imp</i>	Carbapenem
	<i>Escherichia coli</i>	<i>kpc</i>	Carbapenem
	<i>Enterobacter cloacae</i> complex	<i>mecA</i>	Oxacillin
	<i>Enterobacter aerogenes</i>	<i>mecC</i>	Oxacillin
	<i>Proteus</i> spp.	<i>mecC</i>	Oxacillin
	<i>Klebsiella pneumoniae</i>	<i>ndm</i>	Carbapenem
	<i>Klebsiella oxytoca</i>	<i>oxa-23</i>	Carbapenem
	<i>Klebsiella varicola</i>	<i>oxa-24</i>	Carbapenem
	<i>Serratia marcescens</i>	<i>oxa-48</i>	Carbapenem
	<i>Morganella morganii</i>	<i>oxa-58</i>	Carbapenem
Non-fermenting bacteria	<i>Moraxella catarrhalis</i>	<i>shv</i>	Penicillin
	<i>Pseudomonas aeruginosa</i>	<i>sul1</i>	Sulfonamide
	<i>Acinetobacter baumannii</i> complex	<i>tem</i>	Penicillin
	<i>Stenotrophomonas maltophilia</i>	<i>vim</i>	Carbapenem
Others / Fungi	<i>Pneumocystis jirovecii</i>	<i>gyrA83</i>	Fluoroquinolone
	<i>Haemophilus influenzae</i>	<i>gyrA87</i>	Fluoroquinolone
	<i>Mycoplasma pneumoniae</i>		

## Blood culture panel (BCU)

GROUP	PATHOGEN	GENE	RESISTANCE AGAINST
Gram-positive bacteria	Universal Bacteria	<i>aac(6')aph(2'')</i>	Aminoglycoside
	<i>Staphylococcus aureus</i>	<i>ermA</i>	Macrolide/Lincosamide
	Coagulase negative staphylococci	<i>mecA</i>	Oxacillin
	<i>Streptococcus</i> spp.	<i>mecC</i> (LGA251)	Oxacillin
	<i>Streptococcus agalactiae</i>	<i>vanA</i>	Vancomycin
	<i>Streptococcus pneumoniae</i>	<i>vanB</i>	Vancomycin
	<i>Streptococcus pyogenes</i> / <i>dysgalactiae</i>		
	<i>Enterococcus</i> spp.		
	<i>Enterococcus faecalis</i>		
	<i>Listeria monocytogenes</i>		
Enterobacteriaceae	Corynebacteriaceae	<i>aacA4</i>	Aminoglycoside
	<i>Citrobacter freundii</i> / <i>koseri</i>	<i>ctx-M</i>	3rd generation Cephalosporine
	<i>Escherichia coli</i>	<i>kpc</i>	Carbapenem
	<i>Enterobacter cloacae</i> complex	<i>imp</i>	Carbapenem
	<i>Enterobacter aerogenes</i>	<i>ndm</i>	Carbapenem
	<i>Klebsiella oxytoca</i>	<i>oxa-23</i>	Carbapenem
	<i>Klebsiella pneumoniae</i>	<i>oxa-24/40</i>	Carbapenem
	<i>Klebsiella varicola</i>	<i>oxa-48</i>	Carbapenem
	<i>Proteus</i> spp.	<i>oxa-58</i>	Carbapenem
	<i>Serratia marcescens</i>	<i>vim</i>	Carbapenem
Non-fermenting bacteria	<i>Acinetobacter baumannii</i> complex		
	<i>Pseudomonas aeruginosa</i>		
	<i>Stenotrophomonas maltophilia</i>		
other Gram-negative bacteria	<i>Haemophilus influenzae</i>		
	<i>Neisseria meningitidis</i>		
Anaerobic bacteria	<i>Propionibacterium acnes</i>		

\* *Mycobacterium* spp.; Fungi




# 'Sample-in, answer-out'? Evaluation and comprehensive analysis of the Unyvero P50 pneumonia assay

Personne Y. *et al.*, DMID; 86: 2017

London, UK (two hospitals)

Dec 2014 - June 2015

## Pneumonia (42% in ICU)



Sens/Spec: 96% and 33%

## Detection of ARGs

### Routine

36% of strains fully susceptible  
60% resistant to  $\geq 1$  antibiotic  
40% were MDR

### Unyvero

71% of samples with  $\geq 1$  ARG  
(e.g., *bla*<sub>TEM</sub>, *ermB*, *sul1*)

Highly prevalent in commensals

## Species ID

Target Organism	Routine laboratory	Unyvero P50	True Positive (Routine and Unyvero P50)	False Positive (Unyvero P50 only)	False Negative (Routine only)
<i>A. baumannii</i>	3	10	3	7	0
<i>Enterobacter spp</i>	3	9	3	6	0
<i>E. coli</i>	5	8	5	3	0
<i>H. influenzae</i>	3	7	1	6	2
<i>K. pneumoniae</i>	3	11	3	8	0
<i>M. catarrhalis</i>	3	6	3	3	0
<i>M. morgani</i>	0	1	0	1	0
<i>Proteus spp</i>	1	5	1	4	0
<i>P. aeruginosa</i>	13	19	13	6	0
<i>S. marcescens</i>	6	9	6	3	0
<i>S. aureus</i>	5	11	5	6	0
<i>S. maltophilia</i>	6	27	6	21	0
<i>S. mitis group*</i>	0	13	0	13	0
<i>L. pneumophila</i>	0	0	0	0	0
<i>C. pneumoniae</i>	0	0	0	0	0
<i>P. jirovecii</i>	0	0	0	0	0
<i>K. oxytoca</i>	0	0	0	0	0
<i>E. faecalis</i>	1	N/A	0	0	1
<i>C. koseri</i>	1	N/A	0	0	1
Negative specimens	42	16	N/A	N/A	N/A

Avg. pathogens 0.59 vs. 1.59

### Routine

5.5% polymicrobial



### Unyvero

49% polymicrobial




It is not a quantitative test

	ESBL producer	MRSA	Fluoroquinolone resistance	Carbapenemase producer	AmpC producer
<b>Routine Microbiology</b>	not detected	n = 1	n = 3 1x <i>P. aeruginosa</i> (GyrA 83), 2x <i>E. coli</i> (GyrA 83; GyrA 83 + GyrA 87)	n = 4 3 <i>A. baumannii</i> ( <i>bla</i> <sub>OXA-23</sub> ) 1 <i>P. aeruginosa</i> (no enzyme found)	n = 5 3x <i>S. marcescens</i> 2x <i>E. aerogenes</i> Presumed chromosomal AmpC upregulation
<b>Unyvero P50</b>	n = 1 <i>bla</i> <sub>CTX-M</sub>	n = 3*	n = 3 1x <i>P. aeruginosa</i> (GyrA83, ParC) 2x <i>E. coli</i> (GyrA83, GyrA83 + GyrA87)	n = 5 <i>bla</i> <sub>OXA-51</sub> 2x <i>A. baumannii</i> 1x <i>A. baumannii</i> + <i>S. maltophilia</i> 2x <i>A. baumannii</i> + <i>S. maltophilia</i> + <i>S. aureus</i>	n = 5 2x <i>bla</i> <sub>OXA</sub> 1x <i>M. morgani</i> + <i>S. marcescens</i> 1x <i>P. aeruginosa</i> + <i>S. maltophilia</i> 3x <i>bla</i> <sub>EB</sub> 2x <i>Enterobacter spp.</i> 1x <i>Enterobacter spp.</i> + <i>M. catarrhalis</i>
<b>Concordance</b>	No	1/3	2/3	No	No

# Abbott IRIDICA PCR-ESI MS system



## Ibis T5000, Abbott PLEX-ID

780 Bacteria and Candida, 4 Resistance Markers	Whole blood	IRIDICA BAC BSI Assay	
	Sterile fluids and tissues	IRIDICA BAC SFT Assay	
	BAL*, ETA**	IRIDICA BAC LRT Assay	

KPC, VanA, VanB, mecA

# The IRIDICA BAC BSI Assay: Rapid, Sensitive and Culture-Independent Identification of Bacteria and *Candida* in Blood

Metzgar D. *et al.*,  
PlosOne; July 2016

Johns Hopkins Hospital; Baltimore, USA  
285 whole blood samples (with SIRS)  
Samples processed by Ibis Bioscience

## Limit of detection (LOD) in blood:

*S. aureus*, 32 CFU/mL; *E. faecium* 16 CFU/mL; *K. pneumoniae* 32 CFU/mL

273 (95.7%) valid results

61 negative controls: all negative by IRIDICA

Organisms detected: IRIDICA, n=85; culture, n=45

11 of these 46 pathogens gave subsequent infection

### Species ID

	Matched Positive	BAC BSI Assay + / Culture –	BAC BSI Assay – / Culture +	Matched Negative <sup>B</sup>
Gram-positive (including Mycoplasma)	15	11 <sup>A(5)</sup>	2	207
Gram-negative	13	21 <sup>A(4)</sup>	3	207
Unidentified bacteria	0	1	0	207
Yeast	2	1	0	207
Potential Contaminants (details in Table 4)	2	12 <sup>A(2)</sup>	3	207
Other reportable organisms excluding potential contaminants (n = 550)	0	0	0	207

32 out of 40 by culture (80% agreement)  
Excluding contaminants 30/35 (86% agreement)

### Antibiotic R markers

	Matched Positive	BAC BSI Assay + / Culture –	BAC BSI Assay – / Culture +	Matched Negative
<i>bla</i> <sub>KPC</sub> (carbapenem resistance <sup>A</sup> )	0	0	0	0
<i>vanA/vanB</i> (vancomycin resistance <sup>A</sup> )	0	0	0	4
<i>mecA</i> (methicillin resistance <sup>A</sup> )	6	0	0	3


# Broad-Range Detection of Microorganisms Directly from Bronchoalveolar Lavage Specimens by PCR/Electrospray Ionization-Mass Spectrometry Ullberg M. *et al.*, PlosOne; Jan 2017

Karolinska University Hospital, Sweden

121 BAL samples

May 2014 - March 2015

Routine: 104+/121 (86%)  
IRIDICA: 118+/121 (98%)



## Pathogenic species

	IRIDICA-positive	
	Total	Confirmed by culture and/or PCR
<b>Gram-positive bacteria</b>		
<i>Staphylococcus aureus</i>	33	27 <sup>a</sup>
<i>Streptococcus pneumoniae</i>	17	6
<i>Corynebacterium pseudodiphtheriticum</i>	1	-
<b>Gram-negative bacteria</b>		
<i>Haemophilus influenzae</i>	20	16
<i>Klebsiella pneumoniae</i>	4	4
<i>Enterobacter cloacae</i> -complex	5 <sup>b</sup>	4
<i>Escherichia coli</i>	2 <sup>c</sup>	2
<i>Escherichia vulneris</i>	1	-
<i>Proteus mirabilis</i>	2	1
<i>Proteus vulgaris</i>	1	-
<i>Klebsiella oxytoca</i>	1	-
<i>Pseudomonas aeruginosa</i>	8	7
<i>Stenotrophomonas maltophilia</i>	2	1
<i>Chryseobacterium indologenes</i>	1	-
<b>Atypical/ fastidious species<sup>d</sup></b>		
<i>Nocardia</i> species	1	1
<i>Bordetella pertussis</i>	1	1
<i>Legionella pneumophila</i>	1	1
<i>Mycoplasma pneumoniae</i>	1	1

15 by routine culture  
2 by routine culture

Not detected  
by routine culture  
(but other methods)


## Detection of resistance determinants

- *mecA* in 21/121 (17%), but only 8/21 were also positive for *S. aureus*
- MRSA not detected by culture and only 1 *S. aureus* had growth
- No BAL with *vanA/B* or *bla<sub>KPC</sub>*

Mostly  
CoNS+

# Identification of bacterial pathogens and antimicrobial resistance directly from clinical urines by nanopore-based metagenomic sequencing Schmidt K. *et al.*, JAC; 72: 2017

Norfolk/Norwich Hospitals, UK  
10 urines (>10<sup>7</sup> CFU/mL) and 5 spiked (MDR-*Ec*)  
**MinION**: rapid, low capital cost, small



## Species ID

	CU5	CU6	Urine spiked with <i>E. coli</i> from CU 6	CU7	CU8	CU9	CU10	Urine spiked with <i>E. coli</i> H141480453 run 1	Urine spiked with <i>E. coli</i> H141480453 run 2	Urine spiked with <i>E. coli</i> H141480453 run 3
Best species match to MinION sequence data	<i>K. pneumoniae</i> CG43	<i>E. coli</i> JJ1886	<i>E. coli</i> JJ1886	<i>E. coli</i> PMV-1	<i>E. coli</i> 536	<i>E. cloacae</i> NCTC 9394	<i>K. pneumoniae</i> CG43	<i>E. coli</i> APEC O78	<i>E. coli</i> K-12	<i>E. coli</i> APEC O78
Best species match to Illumina sequence data	<i>K. pneumoniae</i> MGH 78578	<i>E. coli</i> JJ1886	<i>E. coli</i> JJ1886	<i>E. coli</i> IHE3034	<i>E. coli</i> 536	<i>E. cloacae</i> NCTC 9394	<i>K. pneumoniae</i> CG43	<i>E. coli</i> ST410	<i>E. coli</i> ST410	<i>E. coli</i> ST410

In all cases, MinION correctly identified the pathogen (WIMP in 15 min)

## ARGs (vs. MICs/EUCAST and Illumina HiSeq)

Urine and species	Method*	AMP	AMC	TZP	CTX	CTX/CLO	CTX/CLA	CAZ	CAZ/CLA	FEP	FEP/CLA	FOX	AZT	CIP	AMK	TOB	GEN	STR	TMP
CU9 <i>E. cloacae</i>	MICs	>64	64	4	2	≤0.12	2	1	1	≤0.12	0.12	>64	0.25	≤0.12	1	0.5	0.5	S	S
	MinION	<i>ampC</i> ( <i>bla</i> <sub>CMY</sub> <sup>revs</sup> , <i>bla</i> <sub>ACT48</sub> , 24)																	
	Illumina	<i>ampC</i> ( <i>bla</i> <sub>ACT24</sub> )																	
CU10 <i>K. pneumoniae</i>	MICs	>64	32	>64	>256	256	0.125	128	1	64	≤0.06	16	>64	>8	8	>32	>32	R	R
	MinION	<i>bla</i> <sub>TEM1</sub> (mv), <i>bla</i> <sub>OXA1</sub>							<i>bla</i> <sub>CTX-M15</sub> , <i>bla</i> <sub>SHV28</sub> , <i>bla</i> <sub>LEN42</sub>							<i>aac</i> (6')-Ib-cr, <i>qnrB</i>			
	Illumina	<i>bla</i> <sub>TEM1</sub> , <i>bla</i> <sub>OXA1</sub>							<i>bla</i> <sub>CTX-M15</sub> , <i>bla</i> <sub>SHV28</sub> , <i>bla</i> <sub>LEN42</sub>							<i>aac</i> (6')-Ib-cr, <i>aacA4</i> , <i>aacC2</i> , <i>dfrA14</i>			
																<i>gyrA</i> (83:S-I), <i>parC</i> (80:S-I), <i>aac</i> (6')-Ib-cr, <i>qnrB</i>			

Illumina detected 55 ARGs from cultured strains. MinION detected 51 of them.

But, **allelic variants poorly distinguished** (e.g., mutations in *gyrA*/*parC*; cAmpCs and pAmpCs)

# Conclusions

- |                                     |   |  |
|-------------------------------------|---|--|
| <b>Species identification</b>       | → | <ul style="list-style-type: none"><li>- Some systems do not provide species ID<br/>For rectal swabs ID probably not necessary</li><li>- Only some report many species [Verigene, FilmArray, Unyvero, Iridica, MinION]</li></ul>  |
| <b>Resistance mechanisms</b>        | → | <ul style="list-style-type: none"><li>- <u>Most systems have a limited list of ARGs (e.g., only KPC)</u></li><li>- Only some report many ARGs [Verigene, Unyvero, MinION]</li><li>- <u>Inability to detect new/emerging ARGs</u></li><li>- <u>Inability to recognize variants (e.g., cAmpCs vs. pAmpCs)</u></li><li>- <u>Is the detected ARG expressed?</u></li><li>- <u>Only two link the ARG with the species ID</u> [Xpert, MinION]</li></ul> |
| <b>Kit and automated</b>            | → | <ul style="list-style-type: none"><li>- Almost all are relatively easy to use (from samples to results)</li></ul>  |
| <b>High sensitivity/specificity</b> | → | <ul style="list-style-type: none"><li>- In general they have good sensitivity (but LOD &lt; than culture)</li><li>- <u>Specificity is low for non-sterile specimens (e.g., from LRT)</u><br/>Difficult interpretation (what will we tell to clinicians?)</li></ul>   |
| <b>Clinical samples</b>             | → | <ul style="list-style-type: none"><li>- Most systems are used for rectal swabs or positive BCs</li><li>- Only IRIDICA for whole blood</li><li>- <u>"Not ready" for direct implementation in respiratory samples</u></li></ul>  |
| <b>Rapid (TAT)</b>                  | → | <ul style="list-style-type: none"><li>- From &lt;1 h [colorim., Xpert, FilmArray, Eazyplex] to 4-8 h [Unyvero, Iridica, MinION]</li><li>- <u>The overall outcome is affected in positive (but lack of studies)</u></li></ul>   |
| <b>Cost-effective</b>               | → | <ul style="list-style-type: none"><li>- From ~10 Euro [e.g., colorimetric] to ~600 Euro [e.g., MinION]</li><li>- <u>Who pay?</u> Still difficult to convince institutions</li></ul>  |



# THANK YOU!

## Institute for Infectious Diseases - University of Bern, Switzerland

- Dr. Valentina Donà, PhD (PostDoc)
- Dr. Baharak Babouee-Flury (PostDoc)
- Odette J. Bernasconi (PhD student)
- João Pires (PhD student)
- Regula Tinguely (Lab Tec)
- Sara Kasraian (Lab Tech)



SWISS NATIONAL SCIENCE FOUNDATION

Grant No. 153377 to Andrea Endimiani

\* Hans Sigrist Foundation (Bern, Switzerland)