



Bern - November 22nd, 2017
(13:00-13:45)



Symposium for Infectious Diseases

How to stay ahead in the race against resistant pathogens

Emerging clinical needs not covered with the current automation

Prof. Andrea Endimiani, MD, PhD

Institute for Infectious Diseases (IfIK)

University of Bern, Switzerland


TREATMENT SETTING

- **Rapid species identification and antibiogram**
 - for serious infections (e.g., BSIs)
 - for patients at risk (e.g., traveled/admitted from ≠ countries, in ICU)
 - enable effective antibiotic stewardship interventions
- **Rapid detection of key resistance mechanisms or genes**
 - ESBLs and carbapenemases
 - affecting other key antibiotics (e.g., polymyxins)
- **Rapid MIC values for key antibiotics**
 - Cefepime, carbapenems, polymyxins


HOSPITAL HYGIENE SETTING

- **Rapid screening to identify patients colonized with MDR bacteria**
 - for isolation (prevention of outbreaks, save costs)
- **Rapid identification of clonal MDR strains**
 - identifying sources, tracking dynamics, intervention
- **Microbiological lab data and algorithms ***

Standard culture approach



Even more days to get the final report including resistance mechanism(s)/gene(s)



Empirical treatment

Isolation



Directed treatment

Septic shock

Inappropriate antibiotic use within the first 6 hours is associated with a 5-fold higher mortality risk

[Kumar A et al., Chest, 2009]

BSIs



Mean decrease in survival of 7.6% for each hour after the onset of BSI until initiation of effective antibiotics

[Kumar A et al., Crit Care Med, 2006]

MDR Gram-neg BSI

It costs 30.000 \$ more (11 vs. 6 days of hospitalization)

[Schwaber MJ et al., JAC, 2007]



→ ID/AST
Still
standard

Advantages

Activity/productivity (increase processing of diagnostic samples).

- Improvement of laboratory workflow (dashboards).
- Management reports.
- Cost savings.

Quality and reproducibility.

- Inoculation: Improved yield of isolated colonies.
- Incubation: Improved bacterial growth.

Reduced time to results (ID and AST).

- Decrease hospitalization time, decrease risks of nosocomial infections, treatment improvements.
- Cost savings.

Traceability (barcodes).

- Decrease errors (e.g. sample, media plates, broth switching).

Security.

- Decrease plate transportation.

Labor saving.

- Decrease fastidious and repetitive tasks (e.g. inoculation, plates incubation).
- Release expert staff for added value tasks (e.g. pre- and postanalytic phase, reading, interpretation, troubleshooting, R&D, microscopy).
- Reduce overtime payments.

diagnostic laboratory. Thus, the real benefits of lab automation remain to be demonstrated in objective, comparative and prospective clinical studies performed by independent laboratories and published in peer-reviewed journals. Several

Total Laboratory Automation and MALDI-TOF Improve Turnaround Times in the Clinical

Microbiology Laboratory: a Retrospective Analysis

Theparee T et al., JCM Nov 8, 2017

Bruker MALDI-TOF MS (Oct, 2013) - BD Kiestra TLA (Dec, 2014)

STUDY AIM: TAT for +/- urine cultures before and after implementation

Comparison I:




Impact of Kiestra TLA

| | Time to ID or preliminary negative | | Time to AST or final negative | |
|---------------------------------|------------------------------------|-----------------------------|-------------------------------|-----------------------------|
| | Pre-TLA | Post-TLA | Pre-TLA | Post-TLA |
| Negative cultures (N=42,259) | 17.73 h (14.97-22.25) | 13.62 h (12.60-16.80)*** | 37.38 h (34.35-42.17) † | 38.62 h (36.85-42.53) † |
| Positive cultures (N=17,251) | 18.53 h (15.00-31.62) | 16.92 h (14.95-25.87)*** | 41.80 h (38.08-55.78) | 40.85 h (38.53-48.68)*** |

Comparison II:

Impact of MALDI and/or Kiestra TLA

| | CONV | MALDI | TLA1 | TLA2 |
|-------------|--------------------------|--------------------------|--------------------------|--------------------------|
| Overall | | | | |
| Time to ID | 21.33 h (14.96-36.79) | 18.02 h (14.70-30.70) | 16.45 h (14.67-24.28) | 16.79 h (14.83-23.55) |
| Time to AST | 42.74 h (39.55-57.43) | 42.33 h (39.30-55.38) | 40.70 h (38.65-47.73) | 40.33 h (38.07-47.08) |




Future developments on lab automation

MALDI-TOF MS



SMART INCUBATORS




Automated colony picking for ID and AST (e.g., MaldiTofA and Colibri)



Automated disk diffusion AST preparation

Automated inhibition zone measurements

Expert systems for dynamic AST interpretation




Automated direct screening for resistance of Gram-negative blood cultures using the BD Kiestra WorkCell

C. S. Heather¹ • M. Maley¹ Eur J Clin Microbiol Infect Dis
02 October 2017

BACTEC 187 +BCs for G- → 2 drops in Copan ESwab → Kiestra Inoqula MH plates → Pre-loaded disk stamp → Kiestra Inoqula MH plates → Kiestra ReadA incubators

15 min

Digital image of automated direct screening for resistance (5 hrs)



BD ReadA
Browser software
(Inhib. zone marked manually)

AST *Enterobacteriaceae* at 5 hrs vs. VITEK2

| | Cut-off | Sensitivity (95% CI) | Specificity (95% CI) | Non-susceptible isolates (%) |
|-------------------------|---------|----------------------|----------------------|------------------------------|
| Ampicillin | 13 mm | 98.3 (93.9–99.8) | 98.2 (90.3–99.9) | 67.6 |
| Ceftriaxone | 19 mm | 100 (92.6–100) | 99.2 (95.5–99.9) | 28.2 |
| Piperacillin–tazobactam | 20 mm | 87.0 (66.4–97.2) | 84.7 (77.8–90.2) | 13.8 |
| Meropenem | 21 mm | 93.8 (69.8–99.8) | 94.8 (90.0–97.7) | 9.4 |
| Ciprofloxacin | 17 mm | 96.6 (82.2–99.9) | 95.7 (91.0–98.4) | 17.1 |
| Gentamicin | 14 mm | 97.4 (86.2–99.9) | 100 (97.1–100) | 23.5 |

AST *P. aeruginosa* at 8 hrs vs. VITEK2

| Antibiotic (zone size cut-off) | Sensitivity (95% CI) | Specificity (95% CI) | Non-susceptible isolates (%) |
|---------------------------------|----------------------|----------------------|------------------------------|
| Ampicillin (13 mm) | 100 (73.5–100) | N/A ^a | 100 |
| Ceftriaxone (19 mm) | 91.7 (61.5–99.8) | N/A ^a | 100 |
| Piperacillin–tazobactam (20 mm) | 100 (2.5–100) | 100 (71.5–100) | 7.7 |
| Meropenem (21 mm) | 100 (15.8–100) | 80 (44.4–97.5) | 7.7 |
| Ciprofloxacin (17 mm) | 100 (2.5–100) | 100 (71.5–100) | 7.7 |
| Gentamicin (14 mm) | 100 (2.5–100) | 90.9 (58.7–99.8) | 7.7 |

82.8% were evaluable at 4 hrs

Ideal system

- Species identification (ID)
- AST (MIC) and interpretation (S, I, R)
- Genes and/or resistance mechanisms
- Kit and automated (easy to use)
- High sensitivity/specificity
- Clinical samples (or at least +BCs)
- Rapid
- Cost-effective

Amplex Eazyplex




- Preparation (**5 min**)
- No DNA extraction/purification
- Run time (**<30 min**)

LAMP

Loop-mediated isothermal Amplification

Real-time fluorescent measurement



~55 Euro

eazyplex® MRSA

S. aureus
mecA
mecC
S. epidermidis

eazyplex® VRE

eazyplex® SuperBug mcr-1

eazyplex® SuperBug CRE

KPC,
NDM,
OXA-48
and OXA 181,
VIM, as well as
CTX-M-1 and
CTX-M-9 group

eazyplex® SuperBug complete

Version A

KPC
NDM
OXA-48
VIM
OXA-23 group
OXA-40 group
OXA-58 group

Version B

KPC
NDM
OXA-48
VIM
OXA-23 group
OXA-40 group
OXA-181

Evaluation of loop-mediated isothermal amplification for the rapid identification of bacteria and resistance determinants in positive blood cultures

Rödel J. *et al.*, Eur J Clin Microbiol Infect Dis; 2017

Germany, 2015

370 +BCs (of which 140 G+)

Test preparation from 25 µL of +BC: 5-10 min

eazyplex® MRSA

Table 1 Analytes detected by the eazyplex® assays used in this study^a

| MRSA | GP | GN |
|-----------------------|---------------------------|-------------------------------------|
| <i>S. aureus</i> | <i>E. faecalis</i> | <i>E. coli</i> |
| <i>S. epidermidis</i> | <i>Enterococcus</i> spp. | <i>K. pneumoniae</i> |
| <i>mecA</i> | <i>S. pneumoniae</i> | Enterobacteriaceae |
| <i>mecC</i> | <i>Streptococcus</i> spp. | <i>Pseudomonas</i> spp. |
| | <i>vanA</i> | <i>bla</i> _{CTX-M-1} group |
| | <i>vanB</i> | <i>bla</i> _{CTX-M-9} group |

In-house LAMP assays (not available)

| ID § Resistance | True-positive (n) | True-negative (n) | False-positive (n) | False-negative (n) | No. positive/no. tested (%) (95% CI ^a) | |
|--|----------------------|----------------------|-----------------------|-----------------------|--|------------------|
| | | | | | Sensitivity | Specificity |
| MRSA | | | | | | |
| <i>S. aureus</i> | 31 | 106 | 2 | 0 | 100 (88.8–100) | 98.2 (93.5–99.8) |
| MRSA ^d | 6 | 25 | 0 | 0 | 100 (54.1–100) | 100 (86.3–100) |
| <i>S. epidermidis</i> | 72 | 60 | 1 | 6 | 92.3 (84–97.1) | 98.4 (91.2–100) |
| <i>mecA</i> (CoNS) | 73 | 30 | 0 | 3 | 96 (88.9–99.2) | 100 (88.4–100) |
| GP | | | | | | |
| <i>E. faecalis</i> | 12 | 58 | 1 | 0 | 100 (73.5–100) | 98.3 (90.9–100) |
| <i>Enterococcus</i> spp. | 41 | 27 | 2 | 1 | 97.6 (87.4–99.9) | 93.1 (77.2–99.1) |
| VRE ^e | 3 | 39 | 0 | 0 | 100 (29.2–100) | 100 (91–100) |
| <i>S. pneumoniae</i> | 2 | 69 | 0 | 0 | 100 (15.8–100) | 100 (94.8–100) |
| <i>Streptococcus</i> spp. | 20 | 50 | 0 | 1 | 95.2 (76.2–99.9) | 100 (92.9–100) |
| GN | | | | | | |
| <i>E. coli</i> | 97 | 58 | 2 | 0 | 100 (96.3–100) | 96.7 (88.5–99.6) |
| <i>E. coli</i> CTX-M ^f | 17 | 80 | 0 | 0 | 100 (89.5–100) | 100 (95.5–100) |
| <i>E. coli</i> ESBL ^g | 17 | 78 | 0 | 2 | 89.5 (66.9–98.7) | 100 (95.3–100) |
| <i>K. pneumoniae</i> | 11 | 143 | 2 | 1 | 91.7 (61.5–99.8) | 98.6 (95.1–99.8) |
| <i>K. pneumoniae</i> ESBL ^h | 4 | 7 | 0 | 0 | 100 (39.8–100) | 100 (59–100) |
| Enterobacteriaceae | 134 | 12 | 11 | 2 | 98.5 (94.8–99.8) | 52.2 (30.6–73.2) |
| <i>Pseudomonas</i> spp. | 9 | 148 | 0 | 0 | 100 (66.4–100) | 100 (97.5–100) |


| Time to results | Threshold time [min; mean values (SD)] | |
|--|--|----------------|
| | True-positive | False-positive |
| MRSA | | |
| <i>S. aureus</i> | 7 (1.5) | 17 |
| <i>mecA</i> (<i>S. aureus</i>) | 11.5 (4) | – |
| <i>S. epidermidis</i> | 11.75 (2.75) | 18.75 |
| <i>mecA</i> (CoNS) | 9.75 (2.5) | – |
| IC | 7.5 (1) | – |
| GP | | |
| <i>E. faecalis</i> | 6.75 (1) | 17.5 |
| <i>Enterococcus</i> spp. | 10 (4) | 17.5 |
| <i>vanA</i> | 7 (1) | – |
| <i>vanB</i> | – | – |
| <i>S. pneumoniae</i> | 7 | – |
| <i>Streptococcus</i> spp. | 13.5 (3.25) | – |
| IC | 9.5 (0.75) | – |
| GN | | |
| <i>E. coli</i> | 7.25 (1.75) | 15 |
| <i>K. pneumoniae</i> | 13.75 (2.5) | 10.75 |
| <i>bla</i> _{CTX-M-1} group ^a | 8.5 (4) | – |
| <i>bla</i> _{CTX-M-9} group ^b | 6.25 (1) | – |
| Enterobacteriaceae | 6.25 (1.75) | 15.75 (2.5) |
| <i>Pseudomonas</i> spp. | 8.25 (1.5) | – |
| IC | 9 (1.75) | – |

From +BC ~30 min

Verigene, Nanosphere



- Load cartridge, consumables, and sample (**5 min**)
- Automated sample preparation and processing
- Place slide from cartridge in reader (**2.5 hrs**)



Microarray approach by using Au-nanoprobe as reporter and silver reduction to enhance signal

Gram-negatives cartridge

| Species | Genus | Resistance |
|-------------------------------|---------------------------|---------------------|
| <i>Escherichia coli</i> * | <i>Acinetobacter</i> spp. | CTX-M (ESBL) |
| <i>Klebsiella pneumoniae</i> | <i>Citrobacter</i> spp. | IMP (carbapenemase) |
| <i>Klebsiella oxytoca</i> | <i>Enterobacter</i> spp. | KPC (carbapenemase) |
| <i>Pseudomonas aeruginosa</i> | <i>Proteus</i> spp. | NDM (carbapenemase) |
| <i>Serratia marcescens</i> | | OXA (carbapenemase) |
| | | VIM (carbapenemase) |

Gram-positives cartridge

| Species | Genus | Resistance |
|--------------------------------------|----------------------------|---------------------------|
| <i>Staphylococcus aureus</i> | <i>Staphylococcus</i> spp. | <i>mecA</i> (methicillin) |
| <i>Staphylococcus epidermidis</i> | <i>Streptococcus</i> spp. | <i>vanA</i> (vancomycin) |
| <i>Staphylococcus lugdunensis</i> | <i>Micrococcus</i> spp. | <i>vanB</i> (vancomycin) |
| <i>Streptococcus anginosus</i> Group | <i>Listeria</i> spp. | |
| <i>Streptococcus agalactiae</i> | | |
| <i>Streptococcus pneumoniae</i> | | |
| <i>Streptococcus pyogenes</i> | | |
| <i>Enterococcus faecalis</i> | | |
| <i>Enterococcus faecium</i> | | |

~50 Euro

Performance Evaluation of the Verigene Gram-Positive and Gram-Negative Blood Culture Test for Direct Identification of Bacteria and Their Resistance Determinants from Positive Blood Cultures in Hong Kong

Gilman K.H. Siu *et al.*, Plos One, Oct 2015

Multicenter study (4 hospitals), 2014

364 BCs (114 Gram-pos; 250 Gram-neg)

ID for Gram-positives (agreement 89.6%)

| Organisms | No. (%) of isolates | | | | | | Sensitivity (%) | Specificity (%) |
|----------------------------|---------------------|----------------------|--------------|---------------|---------|--|-----------------|-----------------|
| | Total | Correctly identified | Not detected | Misidentified | No Call | | | |
| <i>Staphylococcus</i> spp. | 73 (60.3) | 69 (94.5) | - | - | 4 | | 100 | 100 |
| <i>S. aureus</i> | 48 (39.7) | 47 (97.9) | - | - | 1 | | 100 | 100 |
| <i>S. epidermidis</i> | 4 (3.3) | 4 (100) | - | - | - | | 100 | 100 |
| <i>Streptococcus</i> spp. | 27 (22.3) | 25 (92.6) | - | - | 2 | | 100 | 98.9 |
| <i>S. pyogenes</i> | 2 (1.7) | 2 (100) | - | - | - | | 100 | 100 |
| <i>S. agalactiae</i> | 3 (2.5) | 3 (100) | - | - | - | | 100 | 100 |
| <i>S. anginosus</i> gr. | 1 (0.8) | 1 (100) | - | - | - | | 100 | 100 |
| <i>S. pneumoniae</i> | 7 (5.8) | 7 (100) | - | - | - | | 100 | 99.1 |
| <i>E. faecalis</i> | 7 (5.8) | 4 (57.1) | 2 (28.6) | - | 1 | | 66.7 | 100 |
| <i>E. faecium</i> | 5 (4.1) | 2 (40) | 2 (40) | - | 1 | | 50 | 100 |
| <i>Listeria</i> | 1 (0.8) | 1 (100) | - | - | - | | 100 | 100 |

ID for Gram-negatives (agreement 90.5%)

| Organisms | No. (%) of isolates | | | | | | Sensitivity (%) | Specificity (%) |
|---------------------------|---------------------|----------------------|----------------------|----------------|---------|--|-------------------|-----------------|
| | Total | Correctly identified | Not detected | Misidentified | No Call | | | |
| <i>E. coli</i> | 165 (59.1) | 158 (95.8) | 6 (3.6) | 1 ^c | - | | 95.8 | 100 |
| <i>K. pneumoniae</i> | 52 (18.6) | 36 (69.2) | 14 (27) ^d | 2 ^e | - | | 69.2 ^d | 100 |
| <i>P. aeruginosa</i> | 17 (6.1) | 13 (76.5) | 3 (17.6) | - | 1 | | 81.3 | 100 |
| <i>Proteus</i> spp. | 10 (3.6) | 9 (90) | - | - | 1 | | 100 | 100 |
| <i>Enterobacter</i> spp. | 10 (3.6) | 8 (80) | - | - | 2 | | 100 | 99.3 |
| <i>Acinetobacter</i> spp. | 4 (1.4) | 4 (100) | - | - | - | | 100 | 100 |
| <i>K. oxytoca</i> | 3 (1.1) | - | 3 (100) | - | - | | 0 | 99.3 |
| <i>Citrobacter</i> spp. | 1 (0.4) | 1 (100) | - | - | - | | 100 | 100 |

A. hydrophila, *M. morgani*, *Salmonella* spp., *A. faecalis*, *B. pseudomallei*, *H. influenzae*, *R. planticola*, *P. putida*, *S. maltophilia* not detectable with the system (n=16; 5.5%)

MDR organisms

| Drug resistant Organisms | No. of isolates | | | | Sensitivity | Specificity (%) |
|--|-----------------|--------------------|--------------|---------|-------------|-----------------|
| | Total | Correctly Detected | Not Detected | No Call | | |
| Gram Positive | | | | | | |
| MRSA | 27 | 26 | 0 | 1 | 100 | 100 |
| MRSE | 4 | 4 | 0 | 0 | 100 | 100 |
| VRE | 1 | 1 | 0 | 0 | 100 | 100 |
| Gram Negative | | | | | | |
| Cefotaxime resistant Enterobacteriaceae (including ESBL producers) | 61 | 38 | 22 | 1 | 63.3 | 100 |
| ESBL producing Enterobacteriaceae | 46 | 38 | 7 | 1 | 84.4 | 100 |
| MDR Acinetobacter | 3 | 3 | 0 | 0 | 100 | 100 |
| Carbapenem-resistant Pseudomonas | 2 | 0 | 2 | 0 | 0 | 100 |
| Total | 98 | 72 | 24 | 2 | 75.0 | 100 |

Time to Results (40-99 hrs faster than routine)

| Organisms | No. of Isolates ^b | Average Time to identification by Culture-Based Method (h) | Average Time to Result by Verigene Test (h) | Δ Time to Result ^a | | |
|----------------------------|------------------------------|--|---|-------------------------------|--------------|-----------|
| | | | | Average (h) | Range (h) | p-value |
| <i>Staphylococcus</i> spp. | 23 | 63.76 | 2.35 | 61.41 | 22.65-494.65 | p<0.001 |
| MSSA | 7 | 49 | 2.35 | 46.65 | 22.65-79.65 | p = 0.012 |
| MRSA | 10 | 93.45 | 2.35 | 91.1 | 23.65-494.65 | p<0.001 |
| CNS | 8 | 57.19 | 2.35 | 54.84 | 29.65-97.65 | p<0.001 |
| <i>Streptococcus</i> spp. | 13 | 63.86 | 2.35 | 61.51 | 22.65-128.15 | p<0.001 |
| <i>S. pneumoniae</i> | 3 | 42.83 | 2.35 | 40.48 | 22.65-56.65 | p = 0.017 |
| β-haemolytic Strept. | 4 | 71.62 | 2.35 | 69.27 | 33.65-128.15 | p<0.001 |
| Viridians Group | 6 | 55.08 | 2.35 | 52.73 | 25.65-83.65 | p<0.001 |
| <i>Enterococcus</i> spp. | 2 | 78.5 | 2.35 | 76.15 | 53.15-99.15 | p<0.001 |
| VSE | 1 | 55.5 | 2.35 | 53.15 | n/a | n/a |
| VRE | 1 | 101.5 | 2.35 | 99.15 | n/a | n/a |
| Enterobacteriaceae | 79 | 45.98 | 1.88 | 44.1 | 18.12-104.12 | p<0.001 |
| <i>P. aeruginosa</i> | 7 | 52.51 | 1.88 | 50.63 | 28.62-94.62 | p<0.001 |
| <i>Acinetobacter</i> spp. | 1 | 76 | 1.88 | 74.12 | n/a | n/a |

Clinical Impact of Laboratory Implementation of Verigene BC-GN Microarray-Based Assay for Detection of Gram-Negative Bacteria in Positive Blood Cultures

Walker T. *et al.*, J Clin Microbiol; 54:7, 2016

Los Angeles
May 2013 - Nov 2013
Dec 2014 - May 2015

6-month period before Verigene (n=98) vs. 6-month period after Verigene (n=97)

Inpatient stay = 1,400 USD/day

ICU stay = 2,800 USD/day

Test cost (all incl.) = 99 USD

Clinical outcome

| | Pre-BC-GN | Post-BC-GN | P value |
|---|-----------------|------------|---------------------|
| Mean time from initial Gram stain to BC-GN identification, h | NA ^a | 3.5 | NA |
| Mean time from initial Gram stain to organism identification, h | 37.9 | 10.9 | <0.001 ^b |
| Mean time from initial Gram stain to effective therapy, h | | | |
| All cases | 10.2 | 6.5 | 0.12 ^b |
| Cases on suboptimal empirical therapy | 30.3 | 19.1 | 0.12 ^b |
| No. of cases in which therapy was de-escalated | 33 | 37 | 0.66 ^c |
| Mean time from initial Gram stain to de-escalation, h | 40.9 | 34.1 | 0.14 ^b |
| Recurrence of bacteremia, no. (%) | 8 (8.2) | 3 (3.1) | 0.21 ^c |
| Mean total length of stay in hospital, days | 15.2 | 18.0 | 0.52 ^b |
| Mean length of hospital stay after positive culture, days | 9.7 | 9.4 | 0.87 ^b |
| Mean length of stay in ICU, days | 16.2 | 12.0 | 0.03 ^b |
| 30-day mortality, no. (%) | 19 (19.2) | 8 (8.1) | 0.04 ^c |
| ESBL cases, no. | 15 | 11 | 0.53 ^c |
| Length of stay in hospital, days | 12.0 | 13.5 | 0.59 ^b |
| Mean time to effective therapy, h (no.) ^d | 41.4 (9) | 7.3 (9) | 0.04 ^b |
| 30-day mortality, no. (%) | 4 (26.7) | 0 (0) | 0.11 ^c |

Saved 11,661 USD/case ←

Factors associated to length of ICU stay

| Factor | Univariate analysis ^b | | | Multivariate analysis ^b | | |
|------------------------------|----------------------------------|-----------|---------|------------------------------------|-----------|---------|
| | HR | 95% CI | P value | HR | 95% CI | P value |
| Intervention | 1.62 | 1.03–2.54 | 0.04 | 1.79 | 1.07–2.98 | 0.03 |
| Preinfection length of stay | 0.97 | 0.95–0.99 | 0.001 | 0.97 | 0.95–0.98 | 0.0001 |
| Age | 1.01 | 0.99–1.03 | 0.22 | 1.0 | 0.98–1.01 | 0.90 |
| Sex (female) | 0.64 | 0.40–1.03 | 0.06 | 0.69 | 0.41–1.15 | 0.15 |
| Time to effective therapy | 1.01 | 0.78–1.32 | 0.94 | 1.02 | 0.75–1.38 | 0.90 |
| Disease severity (APACHE II) | 1.00 | 0.97–1.02 | 0.74 | 0.97 | 0.93–1.0 | 0.06 |
| Mortality | 1.03 | 0.57–1.65 | 0.91 | 1.43 | 0.83–2.48 | 0.20 |



Multiv. log. regression for association with 30-d mortality

| Factor | Odds ratio | 95% CI | P value |
|---------------------|------------|-----------|---------|
| Intervention | 0.81 | 0.67–0.98 | 0.03 |
| Length of stay, ICU | 1.00 | 1.00–1.00 | 0.94 |
| Age | 1.01 | 1.00–1.01 | 0.11 |
| Sex (female) | 1.07 | 0.88–1.30 | 0.48 |
| APACHE II | 1.01 | 0.99–1.02 | 0.36 |


BioFire FilmArray



- Preparation of the pouch
- Add pouch to FilmArray station (overall, **2 min**)
- Run time of about **1 h**



Magnetic beads



Multiplex PCR
(RT-PCR for RNA target)

Multiplex real-time PCR
and Melting analysis

FilmArray® BCID Panel targets:

~100 Euro

| Gram+ Bacteria | Gram- Bacteria |
|--|--|
| <i>Enterococcus</i> <i>Listeria monocytogenes</i> <i>Staphylococcus</i> <i>Staphylococcus aureus</i> <i>Streptococcus</i> <i>Streptococcus agalactiae</i> <i>Streptococcus pyogenes</i> <i>Streptococcus pneumoniae</i> | <i>Acinetobacter baumannii</i> <i>Haemophilus influenzae</i> <i>Neisseria meningitidis</i> <i>Pseudomonas aeruginosa</i> <i>Enterobacteriaceae</i> <i>Enterobacter cloacae</i> complex <i>Escherichia coli</i> <i>Klebsiella oxytoca</i> <i>Klebsiella pneumoniae</i> <i>Proteus</i> <i>Serratia marcescens</i> |
| Yeast | Antibiotic Resistance |
| <i>Candida albicans</i> <i>Candida glabrata</i> <i>Candida krusei</i> <i>Candida parapsilosis</i> <i>Candida tropicalis</i> | mecA - methicillin resistance vanA/B - vancomycin resistance KPC - carbapenem resistance |

Species ID

| Species ID | Isolates detected ^a : BCID/comparator | | No. of results: BCID/comparator | | | | Sensitivity or PPA ^b : TP/(TP + FN) (%) | Specificity or NPA ^b : TN/(TN + FP) (%) |
|--------------------------------|---|---------------|------------------------------------|-----------|-----------|-----------|---|---|
| | Clinical arm | Seeded arm | TP +/+ | FP +/- | FN -/+ | TN -/- | | |
| Gram-positive bacteria | | | | | | | | |
| <i>Enterococcus</i> | 102/101 | 29/29 | 127 | 4 | 3 | 2,073 | 127/130 (97.7) | 2,073/2,077 (99.8) |
| <i>L. monocytogenes</i> | 0/0 | 36/36 | 36 | 0 | 0 | 2,171 | 36/36 (100) | 2,171/2,171 (100) |
| <i>Staphylococcus</i> | 780/797 | 2/1 | 770 | 12 | 28 | 1,397 | 770/798 (96.5) | 1,397/1,409 (99.1) |
| <i>S. aureus</i> | 257/257 | 0/0 | 253 | 4 | 4 | 1,946 | 253/257 (98.4) | 1,946/1,950 (99.8) |
| <i>Streptococcus</i> | 140/141 | 63/62 | 198 | 5 | 5 | 1,999 | 198/203 (97.5) | 1,999/2,004 (99.8) |
| <i>S. agalactiae</i> (group B) | 18/18 | 18/18 | 36 | 0 | 0 | 2,171 | 36/36 (100) | 2,171/2,171 (100) |
| <i>S. pneumoniae</i> | 26/25 | 12/12 | 36 | 2 | 1 | 2,168 | 36/37 (97.3) | 2,168/2,170 (99.9) |
| <i>S. pyogenes</i> (group A) | 8/7 | 31/31 | 38 | 1 | 0 | 2,168 | 38/38 (100) | 2,168/2,169 (99.9) |
| Total | 1,331/1,346 | 191/189 | 1,494 | 28 | 41 | 16,093 | 1,494/1,535 (97.3) | 16,093/16,121 (99.8) |

Gram-negative bacteria

| | | | | | | | | | |
|---------------------------|---------|---------|------|----|----|--------|---------------------------|----------------------|---------------------------|
| <i>A. baumannii</i> | 16/14 | 40/37 | 51 | 5 | 0 | 2,151 | 51/51 (100) | 2,151/2,156 (99.8) | <i>R. ornithinolytica</i> |
| <i>Enterobacteriaceae</i> | 307/310 | 187/188 | 490 | 4 | 8 | 1,705 | 490/498 (98.4) | 1,705/1,709 (99.8) | |
| <i>E. cloacae</i> complex | 24/22 | 17/17 | 38 | 3 | 1 | 2,165 | 38/39 (97.4) | 2,165/2,168 (99.9) | |
| <i>E. coli</i> | 149/148 | 6/5 | 150 | 5 | 3 | 2,049 | 150/153 (98.0) | 2,049/2,054 (99.8) | |
| <i>K. oxytoca</i> | 6/6 | 54/58 | 59 | 1 | 5 | 2,142 | 59/64 (92.2) ^c | 2,142/2,143 (99.9) | |
| <i>K. pneumoniae</i> | 74/71 | 37/34 | 102 | 9 | 3 | 2,093 | 102/105 (97.1) | 2,093/2,102 (99.6) | |
| <i>Proteus</i> | 22/22 | 17/17 | 39 | 0 | 0 | 2,168 | 39/39 (100) | 2,168/2,168 (100) | |
| <i>S. marcescens</i> | 22/22 | 55/55 | 76 | 1 | 1 | 2,129 | 76/77 (98.7) | 2,129/2,130 (99.9) | |
| <i>H. influenzae</i> | 8/8 | 35/35 | 43 | 0 | 0 | 2,164 | 43/43 (100) | 2,164/2,164 (100) | |
| <i>N. meningitidis</i> | 1/1 | 35/35 | 36 | 0 | 0 | 2,171 | 36/36 (100) | 2,171/2,171 (100) | |
| <i>P. aeruginosa</i> | 52/52 | 0/0 | 51 | 1 | 1 | 2,154 | 51/52 (98.1) | 2,154/2,155 (99.9) | |
| Total | 681/676 | 483/481 | 1135 | 29 | 22 | 23,091 | 1,135/1,157 (98.1) | 23,091/23,120 (99.9) | |

Evaluation of the FilmArray Blood Culture Identification Panel: Results of a Multicenter Controlled Trial

Salimnia H.. *et al.*, J Clin Microbiol; 54:3, 2016

8 centers in USA, 2012-2014

2,207 BC samples

- 1,568 clinical

- 639 seeded

MDR Organisms

| MDR Organisms | Isolates detected: BCID/comparator | | No. of results: BCID/comparator | | | | Sensitivity or PPA ^a : TP/(TP + FN) (%) | Specificity or NPA ^a : TN/(TN + FP) (%) |
|---|---------------------------------------|---------------|------------------------------------|-----------|-----------|-----------|---|---|
| | Clinical arm | Seeded arm | TP +/+ | FP +/- | FN -/+ | TN -/- | | |
| Antimicrobial resistance gene(s) | | | | | | | | |
| <u>mecA</u> in association with: | | | | | | | | |
| All <i>Staphylococcus</i> isolates detected ^b | 491/494 | 2/2 | 488 | 5 | 8 | 281 | 488/496 (98.4) | 281/286 (98.3) |
| <i>Staphylococcus</i> and <i>S. aureus</i> isolates detected | 137/139 | 0/0 | 137 | 0 | 2 | 118 | 137/139 (98.6) | 118/118 (100) |
| <u>vanA/B</u> in association with | | | | | | | | |
| <i>Enterococcus</i> isolates detected | 36/36 | 28/28 | 64 | 0 | 0 | 67 | 64/64 (100) ^c | 67/67 (100) |
| <u>bla_{KPC}</u> in association with | | | | | | | | |
| <i>Enterobacteriaceae</i> and/ or <i>A. baumannii</i> and/ or <i>P. aeruginosa</i> isolates detected | 6/6 ^d | 33/33 | 39 | 0 | 0 | 558 | 39/39 (100) ^e | 558/558 (100) |

Clinical and economic impact of antimicrobial stewardship interventions with the FilmArray blood culture identification panel *Pardo J et al., DMID; 84, 2016*


Gainesville, Florida

Intervention group (FilmArray)
(Aug 2013 - Jan 2014; n=84)

vs.

Historical control group
(Jan 2012 - Jun 2013; n=252)

55% CoNS



Time to discharge after
CoNS-contaminated BCs

| Variable | All Patients (N = 336) | | | Discharged Within 6 Days of CoNS Contamination (N = 100) | | |
|---------------------------------|-------------------------|---------------------|-------|--|---------------------|-------|
| | Control group (N = 252) | BCID group (N = 84) | P | Control group (N = 75) | BCID group (N = 25) | P |
| Median LOS (days [IQR]) | 7.9 (3.8–16.3) | 7.4 (3–11.7) | 0.178 | 2.9 (2.1–4.4) | 2.3 (1.5–3.1) | 0.008 |
| Median ICU Days (IQR) | 0 (0–5) | 0 (0–1) | 0.026 | 0 (0–0) | 0 (0–0) | 0.838 |
| In-hospital mortality (No. [%]) | 37 (15) | 5 (6) | 0.036 | 4 (6) | 0 (0) | 0.239 |
| Median Costs (\$ [IQR]) | | | | | | |
| Total | 15324 (5517–37305) | 12241 (3080–21027) | 0.030 | 3370 (1537–5190) | 1645 (585–2952) | 0.016 |
| Pharmacy | 2119 (598–5503) | 1538 (366–3824) | 0.073 | 429 (157–712) | 214 (75–361) | 0.010 |
| Ward room | 2806 (0–7427) | 3972 (947–8041) | 0.514 | 1852 (891–2778) | 947 (0–1828) | 0.025 |
| ICU room | 1222 (0–10218) | 0 (0–3756) | 0.008 | 0 (0–0) | 0 (0–0) | 0.659 |
| Laboratory | 708 (241–1804) | 656 (254–1445) | 0.692 | 160 (83–253) | 204 (160–321) | 0.020 |




129 USD for the test

Unnecessary vancomycin use

| Variable | MSSA (N = 32) | | | CoNS Contaminants (N = 184) | | |
|------------------------------|------------------------|--------------------|--------|-----------------------------|---------------------|-------|
| | Control Group (N = 24) | BCID Group (N = 8) | P | Control Group (N = 138) | BCID Group (N = 46) | P |
| Vancomycin started (No. [%]) | 24 (100) | 8 (100) | >0.99 | 100 (72) | 32 (70) | 0.705 |
| 1 dose only (No. [%]) | 0/24 (0) | 4/8 (50) | 0.002 | 12/100 (12) | 7/32 (22) | 0.133 |
| Median duration (h [IQR]) | 66.0 (54.1–78.1) | 14.0 (1–20) | <0.001 | 66.9 (37.6–135) | 39.3 (14.7–96.5) | 0.100 |

Abbott IRIDICA PCR-ESI MS system



1-6 Samples


Setup, Reagent addition, sample addition and lysis per sample is 15min

Nucleic purification with setup is 2hrs 10 min




PCR amplification is 2 hrs 20 min

Desalting is 45 min

Mass Spectrometry with setup is 2 hours (for 1 sample it is about 25 min)



Ibis T5000, Abbott PLEX-ID

| | | | |
|---|----------------------------|-----------------------|---|
| 780 Bacteria and Candida, 4 Resistance Markers | Whole blood | IRIDICA BAC BSI Assay |  |
| | Sterile fluids and tissues | IRIDICA BAC SFT Assay |  |
| | BAL*, ETA** | IRIDICA BAC LRT Assay |  |

KPC, VanA, VanB, mecA

~200-300 Euro



The IRIDICA BAC BSI Assay: Rapid, Sensitive and Culture-Independent Identification of Bacteria and Candida in Blood

Metzgar D. *et al.*,
PlosOne; July 2016

Johns Hopkins Hospital; Baltimore, USA
285 whole blood samples (with SIRS)
Samples processed by Ibis Bioscience

Limit of detection (LOD) in blood: 16-32 CFU/mL

273 (95.7%) valid results

61 negative controls: all negative by IRIDICA

Organisms detected: IRIDICA, n=85; culture, n=45

11 of these 46 pathogens
gave subsequent infection

Species ID

| | Matched Positive | BAC BSI Assay + / Culture – | BAC BSI Assay–/ Culture + | Matched Negative ^B |
|---|------------------|-----------------------------|---------------------------|-------------------------------|
| Gram-positive (including Mycoplasma) | 15 | 11 ^{A(5)} | 2 | 207 |
| Gram-negative | 13 | 21 ^{A(4)} | 3 | 207 |
| Unidentified bacteria | 0 | 1 | 0 | 207 |
| Yeast | 2 | 1 | 0 | 207 |
| Potential Contaminants (details in Table 4) | 2 | 12 ^{A(2)} | 3 | 207 |
| Other reportable organisms excluding potential contaminants (n = 550) | 0 | 0 | 0 | 207 |

32 out of 40 by culture (80% agreement)




Antibiotic R markers

| | Matched Positive | BAC BSI Assay + / Culture – | BAC BSI Assay–/ Culture + | Matched Negative |
|--|------------------|-----------------------------|---------------------------|------------------|
| <i>bla_{KPC}</i> (carbapenem resistance ^A) | 0 | 0 | 0 | 0 |
| <i>vanA/vanB</i> (vancomycin resistance ^A) | 0 | 0 | 0 | 4 |
| <i>mecA</i> (methicillin resistance ^A) | 6 | 0 | 0 | 3 |

T2 Magnetic Resonance (T2MR)



- Sample transferred to T2Dx (1 min)
- No extraction/purification
- Time to results (~3 hrs)
- Limit of detection: 1 CFU/mL



T2Bacteria Kit

- **Gram-negatives:** *E. coli*, *K. pneumoniae*, *P. aeruginosa*, *A. baumannii*
- **Gram-positives:** *S. aureus*, *E. faecium*



Accelerate Pheno™ System (BC Kit)




Fully Automated
Entire process with 1 kit
2 min to load and start
Directly from +BC bottle




ID: <1.5 hrs
AST: <7 hrs

~120 Euro

ID: 16 FISH probes



AST: morphokinetic cell analysis



| Gram-positive | Identification | Ampicillin | Ceftaroline | Erythromycin | Daptomycin | Linezolid | Vancomycin | MRSA (Cefoxitin) | MLSD (Erythromycin-Clinamycin) |
|---------------------------|----------------|------------|-------------|--------------|------------|-----------|------------|------------------|--------------------------------|
| <i>S. aureus</i> | • | | • | • | • | • | • | • | |
| <i>S. lugdunensis</i> | • | | | | | | • | • | • |
| CNS spp. | • | | | | • | | • | • | • |
| <i>E. faecalis</i> | • | • | | | • | • | • | | |
| <i>E. faecium</i> | • | • | | | • | • | • | | |
| <i>Streptococcus</i> spp. | • | | | | | | | | |

| Gram-negative | Identification | Ampicillin-Sulbactam | Piperacillin-Tazobactam | Cefepime | Ceftazidime | Ceftiofur | Ertapenem | Meropenem | Aztreonam | Gentamicin | Tobramycin | Ciprofloxacin | Aztreonam |
|--------------------------|----------------|----------------------|-------------------------|----------|-------------|-----------|-----------|-----------|-----------|------------|------------|---------------|-----------|
| <i>E. coli</i> | • | • | • | • | • | • | • | • | • | • | • | • | • |
| <i>Klebsiella</i> spp. | • | • | • | • | • | • | • | • | • | • | • | • | • |
| <i>Enterobacter</i> spp. | • | | • | • | • | • | • | • | • | • | • | • | • |
| <i>Proteus</i> spp. | • | • | • | • | • | • | • | • | • | • | • | • | • |
| <i>Citrobacter</i> spp. | • | | • | • | • | • | • | • | • | • | • | • | • |
| <i>S. marcescens</i> | • | | • | • | • | • | • | • | • | • | • | • | • |
| <i>P. aeruginosa</i> | • | | • | • | • | | • | • | • | • | • | • | |
| <i>A. baumannii</i> | • | | • | | | | | | • | | | | |

Evaluation of the Accelerate Pheno System for Fast Identification and Antimicrobial Susceptibility Testing from Positive Blood Cultures in Bloodstream Infections Caused by Gram-Negative Pathogens *Marschal M et al., JCM, July 2017*

April 2016 - Oct 2016: Southern Germany

115 BSIs due to Gram-negatives (10 poly)

Routine: BACTEC 9240, MALDI, Vitek2, Etest

ACC started on average 8 h after positive signal

Identification

Correct ID:

102/115 (88.7%)

Correct ID for covered species:

102/105 (97.1%)

e.g., B. fragilis, M. morganii, H. alvei not detected

Accelerate Pheno system AST

| Antimicrobial agent | No. of AST results | No. (%) of category agreements | | | No. (%) of minor discrepancies | No. (%) of major discrepancies | No. (%) of very major discrepancies |
|---------------------|--------------------|--------------------------------|------|------------|--------------------------------|--------------------------------|-------------------------------------|
| | | S | R/I | Total | | | |
| SAM | 66 | 31 | 32 | 63 (95.5) | 0 (0) | 3 (8.8) | 0 (0) |
| TZP | 91 | 78 | 6 | 84 (92.3) | 0 (0) | 7 (8.2) | 0 (0) |
| FEP | 90 | 73 | 6/1* | 80 (88.9) | 7 (7.8) | 3 (3.6) | 0 (0) |
| CRO | 85 | 75 | 8 | 83 (97.6) | 0 (0) | 2 (2.6) | 0 (0) |
| ETP | 85 | 85 | 0 | 85 (100) | 0 (0) | 0 (0) | 0 (0) |
| MEM | 94 | 91 | 2 | 93 (98.9) | 1 (1.1) | 0 (0) | 0 (0) |
| AMK | 95 | 91 | 2 | 93 (97.9) | 1 (1.1) | 0 (0) | 1 (33.3) |
| GEN | 94 | 84 | 9 | 93 (98.9) | 0 (0) | 1 (1.2) | 0 (0) |
| TOB | 94 | 82 | 10 | 92 (97.9) | 1 (1.1) | 1 (1.2) | 0 (0) |
| CIP | 94 | 69 | 21 | 90 (95.7) | 4 (4.3) | 0 (0) | 0 (0) |
| CST | 88 | 85 | 0 | 85 (96.6) | 0 (0) | 3 (3.4) | 0 (0) |
| Total | 976 | 844 | 97 | 941 (96.4) | 14 (1.4) | 20 (2.3) | 1 (1.0) |

| | | MIC (μg/ml) ^c | | | | | | | | | | | |
|---------------------------|-------------------|--------------------------|----------|---------|---------|-----------|-----------|----------|----------|----------|-----------|----------|---------------|
| Isolate (n = 13) | Test ^b | SAM | TZP | FEP | CRO | ETP | MEM | AMK | GEN | TOB | CIP | CST | Comment |
| <i>C. freundii</i> (1) | ACC | – | 128 (R) | ≤1 (S) | ≥8 (R) | 0.25 (S) | ≤0.25 (S) | ≤4 (S) | ≤1 (S) | ≤1 (S) | ≤0.25 (S) | ≤0.5 (S) | AmpC positive |
| | Etest | ND | >256 (R) | 1 (S) | >32 (R) | 0.094 (S) | 0.047 (S) | 1.5 (S) | 0.25 (S) | 1 (S) | 0.047 (S) | 1.5 (S) | |
| <i>E. coli</i> (9) | ACC | 32 (R) | ≤4 (S) | 16 (R) | ≥8 (R) | 0.25 (S) | ≤0.25 (S) | ≤4 (S) | ≤1 (S) | ≤1 (S) | ≤0.25 (S) | ≤0.5 (S) | ESBL positive |
| | Etest | 16 (R) | 0.5 (S) | 6 (R) | >32 (R) | 0.012 (S) | 0.016 (S) | 0.75 (S) | 0.25 (S) | 0.5 (S) | 0.094 (S) | 1 (S) | |
| <i>E. coli</i> (10) | ACC | ≥64 (R) | ≤4 (S) | ≥32 (R) | ≥8 (R) | 0.25 (S) | ≤0.25 (S) | ≤4 (S) | 16 (R) | 8 (R) | ≥8 (R) | 2 (S) | ESBL positive |
| | Etest | 64 (R) | 0.75 (S) | 16 (R) | >32 (R) | 0.023 (S) | 0.032 (S) | 1.5 (S) | 0.5 (S) | 1 (S) | >32 (R) | 1 (S) | |
| <i>P. aeruginosa</i> (11) | ACC | – | 128 (R) | 32 (R) | – | – | ≥16 (R) | ≥128 (R) | ≥32 (R) | ≥32 (R) | ≥8 (R) | ≤0.5 (S) | MDR phenotype |
| | Etest | ND | >256 (R) | 96 (R) | ND | ND | >32 (R) | 24 (R) | >256 (R) | >256 (R) | >32 (R) | 2 (S) | |

Species ID: 1.35 h – 27.49 hrs before culture-based ID





AST: 6.65 h – 40.39 hrs earlier than conventional tests

Check-Direct Screening, Check-Points

**Rectal
swab**
(Copan, ESwab)



- Rapid preparation
- Time to results (<3 hrs)



Check-Direct **ESBL**
Screen for BD MAX™

CTX-M-1 group
CTX-M-2 group
CTX-M-9 group
SHV-ESBL

Check-Direct **CPE**
Screen for BD MAX™

KPC
OXA-48-like
VIM
NDM

~30 Euro

Performance of the BD MAX™ instrument with Check-Direct CPE real-time PCR for the detection of carbapenemase genes from rectal swabs, in a setting with endemic dissemination of carbapenemase-producing *Enterobacteriaceae* Antonelli *et al.*, DMID, 2016

Florence, Italy
557 rectal swabs

Limit of detection (LOD)

| Target | LOD BDMix (CFU/ml) | LOD CPMix (CFU/ml) |
|------------------------------|---------------------|---------------------|
| <i>bla</i> _{KPC} | 9×10 ⁵ | 9×10 ³ |
| <i>bla</i> _{VIM} | 4.5×10 ² | 4.5×10 ² |
| <i>bla</i> _{OXA-48} | 8.5×10 ¹ | 8.5×10 ¹ |
| <i>bla</i> _{NDM} | 7.3×10 ⁴ | 7.3×10 ¹ |

Culture vs. C-D direct CPE (2 ≠ master mix)

Summary of the results

| | | Colonization status | |
|---------------|------------|--------------------------|-------------------|
| | | Positive (n. 29) | Negative (n. 528) |
| CHSM and BERM | Positive | 24 (82.8 %) ^a | 0 |
| | Negative | 5 (17.2 %) | 528 (100 %) |
| BDMix | Positive | 27 (93.1 %) | 14 (2.7 %) |
| | Negative | 2 (6.9 %) | 496 (93.9 %) |
| | Unresolved | - | 18 (3.4 %) |
| CPMix | Positive | 29 (100 %) | 15 (2.9 %) |
| | Negative | 0 | 507 (96.0 %) |
| | Unresolved | - | 6 (1.1 %) |

ChromID CARBA SMART (w/wo broth enrichment) →
[5 samples not detected] *

TAT from 18-24 hrs (direct culture)
or 48 hrs (broth enrichment) to 3 hrs

23 KPC
5 VIM
1 OXA-48

Performance of Check-Direct CPE screening for BD MAX

| | Positives for <i>bla</i> _{KPC} | | Positives for <i>bla</i> _{VIM} | | Positives for <i>bla</i> _{OXA-48} | | Positives for <i>bla</i> _{NDM} | |
|-------------|---|-------|---|-------|--|-------|---|-------|
| | BDMix | CPMix | BDMix | CPMix | BDMix | CPMix | BDMix | CPMix |
| sensitivity | 91,3% | 100% | 100% | 100% | 100% | 100% | - | - |
| specificity | 98,6% | 99,6% | 99,3% | 97,6% | 99,3% | 99,6% | 100% | 100% |
| PIR | 65.21 | 250 | 142.86 | 41.67 | - ^a | - | - | - |
| NLR | 0.09 | 0 | 0 | 0 | - | - | - | - |

Cepheid GeneXpert



- Add aliquot to elution, vortex, transfer to port S
- Insert cartridge to station (overall, **1 min**)
- Run time (**<1 h**)

Real-time multiplex PCR

- Smart fluidic system
- Filtering and Sonication (DNA)
- Fluorescent-labeled hybr. probes (6 colors)
- Internal control

Xpert® Carba-R

KPC
NDM
OXA-48-like
VIM
IMP-1

Xpert® vanA

Xpert® MRSA/SA BC

MSSA
MRSA
CoNS

~50 Euro

*Targets: spa, mecA and
junction SCCmec-orfX*

Multisite Evaluation of Cepheid Xpert Carba-R Assay for Detection of Carbapenemase-Producing Organisms in Rectal Swabs JCM, 54:7; 2016

M. Tato,^a P. Ruiz-Garbajosa,^a M. Traczewski,^b A. Dodgson,^c A. McEwan,^c R. Humphries,^d J. Hindler,^d J. Veltman,^e H. Wang,^f R. Cantón^a

4 centers (2 USA, 1 UK, 1 Spain)
July 2013 - Feb 2014

Results obtained in 32-48 min

| Xpert Carba-R assay result | Clinical specimens (n = 383) | Contrived specimens (n = 250) | All specimens (n = 633) |
|---|------------------------------|-------------------------------|-------------------------|
| Positive (single and/or combined targets) | 42 | 107 | 149 23.5% |
| IMP-1 | 0 | 25 | 25 |
| VIM | 2 | 24 | 26 |
| NDM | 2 | 23 | 25 |
| KPC | 13 | 19 | 32 |
| OXA-48 | 20 | 15 | 35 |
| VIM + OXA-48 | 4 | 0 | 4 |
| NDM + KPC | 1 | 0 | 1 |
| IMP-1 + NDM | 0 | 1 | 1 |
| Negative | 341 | 143 | 484 |

Results by individual targets

| Xpert Carba-R assay | Reference method (culture plus sequencing) | | | | | | Total |
|---------------------|--|-----|-----|-----|--------|----------|-------|
| | IMP-1 | VIM | NDM | KPC | OXA-48 | Negative | |
| IMP-1 | 26 | 0 | 0 | 0 | 0 | 0 | 26 |
| VIM | 0 | 29 | 0 | 0 | 0 | 1 | 30 |
| NDM | 0 | 0 | 26 | 0 | 0 | 1 | 27 |
| KPC | 0 | 0 | 0 | 29 | 0 | 4 | 33 |
| OXA-48 | 0 | 0 | 0 | 0 | 38 | 1 | 39 |
| Negative | 1 | 2 | 0 | 1 | 2 | 3,004 | 3,010 |
| Total | 27 | 31 | 26 | 30 | 40 | 3,011 | 3,165 |

6 cases of multiple genes (4 not detected by ref. method)
6 out of 154 specimens not detected by Xpert

Performance vs. reference method

| Xpert Carba-R assay | Reference method (culture plus sequencing) | | |
|---------------------|--|--------------|-----------|
| | No. positive | No. negative | Total No. |
| Positive | 142 | 7 | 149 |
| Negative | 6 | 478 | 484 |
| Total | 148 | 485 | 633 |

Performance for different targets

| Target gene | Sensitivity (% [95% CI]) | Specificity (% [95% CI]) | PPV (%) | NPV (%) |
|-------------|--------------------------|--------------------------|---------|---------|
| IMP-1 | 96.3 (81.0–99.9) | 100 (99.4–100) | 100 | 99.8 |
| VIM | 93.5 (78.6–99.2) | 99.8 (99.1–100) | 96.7 | 99.7 |
| NDM | 100 (86.8–100) | 99.8 (99.1–100) | 96.3 | 100 |
| KPC | 96.7 (82.8–99.9) | 99.3 (98.3–99.8) | 87.9 | 99.8 |
| OXA-48 | 95.0 (83.1–99.4) | 99.8 (99.1–100) | 97.4 | 99.7 |

Sensitivity: 96.6% - Specificity: 98.6%

Comparison of the Next-Generation Xpert MRSA/SA BC Assay and the GeneOhm StaphSR Assay to Routine Culture for Identification of *Staphylococcus aureus* and Methicillin-Resistant *S. aureus* in Positive-Blood-Culture Broths

Buchan BW *et al.*, J Clin Microb, 53:3; 2015

Multicenter study (8 hospitals, USA)

795 BCs (30%, *S. aureus*; 13% MRSA)

Identification of *Staphylococcus aureus*

| Test | Site | Total no. of specimens tested | No. with result ^a : | | | | Sensitivity (% [95% CI]) ^b | Specificity (% [95% CI]) |
|------------------|-------|-------------------------------|--------------------------------|----------------|-----|----------------|---------------------------------------|--------------------------|
| | | | TP | FP | TN | FN | | |
| Xpert MRSA/SA BC | A | 63 | 19 | 0 | 44 | 0 | 100 (82.4–100) | 100 (92.0–100) |
| | B | 91 | 32 | 0 | 59 | 0 | 100 (89.1–100) | 100 (93.9–100) |
| | C | 44 | 17 | 0 | 27 | 0 | 100 (80.5–100) | 100 (87.2–100) |
| | D | 70 | 26 | 0 | 44 | 0 | 100 (86.8–100) | 100 (92.0–100) |
| | E | 130 | 27 | 2 ^c | 100 | 1 ^d | 96.4 (79.7–99.8) | 98.0 (92.4–99.6) |
| | F | 211 | 65 | 1 ^e | 145 | 0 | 100 (94.5–100) | 99.3 (96.2–100) |
| | G | 126 | 31 | 0 | 95 | 0 | 100 (88.8–100) | 100 (96.2–100) |
| | H | 57 | 18 | 0 | 39 | 0 | 100 (81.5–100) | 100 (91.0–100) |
| | Total | 792 | 235 | 3 | 553 | 1 | 99.6 (97.7–99.9) | 99.5 (98.4–99.9) |
| GeneOhm StaphSR | Total | 782 | 234 | 19 | 527 | 2 | 99.2 (97.0–99.9) | 96.5 (94.6–97.9) |


Statistically equivalent

Identification of MRSA


| Test | Site | Total no. of specimens tested | No. with result ^a : | | | | Sensitivity (% [95% CI]) ^b | Specificity (% [95% CI]) |
|------------------|-------|-------------------------------|--------------------------------|----------------|-----|----------------|---------------------------------------|--------------------------|
| | | | TP | FP | TN | FN | | |
| Xpert MRSA/SA BC | A | 63 | 11 | 0 | 52 | 0 | 100 (71.5–100) | 100 (93.2–100) |
| | B | 91 | 16 | 0 | 75 | 0 | 100 (79.4–100) | 100 (95.2–100) |
| | C | 44 | 7 | 0 | 36 | 1 | 87.5 (47.3–99.7) | 100 (90.3–100) |
| | D | 70 | 11 | 1 ^c | 58 | 0 | 100 (71.5–100) | 98.3 (90.9–100) |
| | E | 130 | 10 | 1 ^d | 119 | 0 | 100 (69.2–100) | 99.2 (95.4–100) |
| | F | 211 | 25 | 1 ^c | 185 | 0 | 100 (86.3–100) | 99.5 (97.0–100) |
| | G | 126 | 19 | 0 | 106 | 1 | 95.0 (73.0–99.7) | 100 (96.2–100) |
| | H | 57 | 4 | 0 | 53 | 0 | 100 (81.5–100) | 100 (91.6–100) |
| | Total | 792 | 103 | 3 | 684 | 2 ^e | 98.1 (93.3–99.8) | 99.6 (98.7–99.9) |
| GeneOhm StaphSR | Total | 782 | 99 | 15 | 662 | 6 ⁱ | 94.3 (88.0–97.9) | 97.8 (96.4–98.8) |

Statistically equivalent

DiversiLab System: rep-PCR




6 hrs




Repetitive extragenic palindromic PCR

1 h




~100 Euro



Characterization of *bla*_{KPC}-containing *Klebsiella pneumoniae* isolates detected in different institutions in the Eastern USA

Endimiani A et al., JAC, 2009



22-33 strains in one day

Dedicated software
Online comparison with
database

Evaluation of the DiversiLab System for Detection of Hospital Outbreaks of Infections by Different Bacterial Species[▽]


A. C. Fluit,^{1*} A. M. Terlingen,¹ L. Andriessen,¹ R. Ikawaty,¹ R. van Mansfeld,¹ J. Top,¹
J. W. Cohen Stuart,¹ M. A. Leverstein-van Hall,^{1,2} and C. H. E. Boel¹

JOURNAL OF CLINICAL MICROBIOLOGY, Nov. 2010,

TABLE 1. Characteristics of DL compared with those of the reference test method

| | Species (no. of isolates) and test method | Typeability (%) | Simpson's index of diversity (95% CI) ^a | Adjusted Rand's coefficient | Wallace's coefficient 1 | Wallace's coefficient 2 |
|-------------|--|-----------------|--|-----------------------------|-------------------------|-------------------------|
| OK ← | <i>E. cloacae</i> (28) DL PFGE | 100 96.4 | 0.794 (0.632–0.955) 0.843 (0.703–0.983) | 0.691 | 0.682 | 0.818 |
| OK ← | <i>Klebsiella</i> spp. (23) DL PFGE | 100 95.7 | 0.964 (0.868–0.991) 0.938 (0.894–0.982) | 0.623 | 0.808 | 1.0 |
| Not good* ← | <i>E. coli</i> (38) DL MLVA | 100 100 | 0.972 (0.952–0.992) 0.966 (0.941–0.990) | 0.390 | 0.450 | 0.375 |
| Not good ← | <i>P. aeruginosa</i> (13) DL PFGE | 100 92.3 | 0.962 (0.919–1.0) 0.985 (0.949–1.0) | 0.488 | 1.0 | 0.333 |
| | <i>P. aeruginosa</i> (52) DL MLST | 100 100 | 0.977 (0.958–0.997) 0.928 (0.828–0.989) | 0.212 | 0.182 | 0.333 |
| OK ← | <i>Acinetobacter</i> spp. (26) DL PFGE | 100 100 | 0.982 (0.966–0.997) 0.985 (0.970–0.999) | 0.908 | 0.833 | 1.0 |
| OK ← | <i>Stenotrophomonas</i> spp. (25) DL PFGE | 100 96 | 0.963 (0.913–1.0) 0.975 (0.937–1.0) | 0.697 | 0.6 | 0.857 |
| Not good ← | <i>E. faecium</i> (26) DL MLVA/esp | 100 100 | 0.929 (0.868–0.991) 0.892 (0.822–0.962) | 0.623 | 0.826 | 0.543 |
| Not good ← | MRSA (50) DL Hospital hygiene classification | 100 100 | 0.913 (0.859–0.968) 0.976 (0.960–0.991) | 0.266 | 0.189 | 0.571 |
| | MRSA (50) DL PFGE | 100 100 | 0.913 (0.859–0.968) 0.970 (0.953–0.986) | 0.261 | 0.198 | 0.568 |
| | MRSA (50) DL <i>spa</i> typing | 100 100 | 0.913 (0.859–0.968) 0.908 (0.858–0.957) | 0.303 | 0.377 | 0.354 |
| | MRSA (50) DL MLVA | 100 100 | 0.913 (0.859–0.968) 0.971 (0.955–0.988) | 0.252 | 0.189 | 0.571 |
| | | | | | | |

Oxford Nanopore Sequencing



~500 Euro

~200-300 Euro

**Rapid real-time long read WGS
(Point of care, POC)**

Emergence of *Klebsiella pneumoniae* co-producing NDM-1, OXA-48, CTX-M-15, CMY-16, QnrA and ArmA in Switzerland

IJAA, 2014



Salome N. Seiffert^{a,b,c}, Jonas Marshall^d, Vincent Perreten^b, Alessandra Carattoli^{a,b,e}, Hansjakob Furrer^d, Andrea Endimiani^{a,*}

Our experience: ~3 hrs for WGS

Rapid Nanopore Sequencing of Plasmids and Resistance Gene Detection in Clinical Isolates

J. Clin. Microbiol. 11 October 2017

Jamie K. Lemon, Pavel P. Khil, Karen M. Frank, and John P. Dekker[#]




Covert dissemination of carbapenemase-producing *Klebsiella pneumoniae* (KPC) in a successfully controlled outbreak: long- and short-read whole-genome sequencing demonstrate multiple genetic modes of transmission

J Antimicrob Chemother 2017; 72: 3025–3034

UK, 2013–14

Rectal swabs/clinical samples

MinION and Illumina




Clonal?


Plasmid?


Tn4401?

(b)
Direct
transmission
network




(c)
Indirect
transmission
network





"Integrated automation"



Rapid, accurate, cost-effective,
information for treating physicians

BETTER OUTCOME

Disadvantages

No laboratory adaptation to automation (e.g. staff shifts, training, 24/7)

- Misuse of tools
- Expectations for increased productivity not achieved

Crash of automat (backup needed).

- Good support and maintenance essential.
- Expensive maintenance budget.

Staff turnover (boring and lonely work?).

- Lab automation needs to be a project that includes everybody.
- Aim is not to replace experienced laboratory technicians but to assist them in their daily tasks.

Only eye is used.

- Smelling or other sensing of colony consistency disappears.
- More difficult to identify unusual/new species.

Security.

- Inoculation of sensitive samples (e.g. sputum, blood culture).
- Contamination of specimen processors and incubators (e.g. fungus spores, biosafety class 3 microorganisms).

Loss of microbiologic knowledge.

- Decrease in analytical variability.
- Standardized microbiologic factory (you find what you are looking for).

Clinical Microbiology and Infection, Volume 22 Number 3, March 2016

to new insights into disease pathogenesis. It is true that we will need a considerable amount of expertise to work with these advanced approaches. However, it is also imperative for us to accept these challenges and to enact wise choices in the selection of diagnostics. We are not certain which method will best suit our needs, but, perhaps, "chance will favor a prepared mind".

Robert A. Bonomo and Andrea Endimiani, DMID, 2013

THANK YOU!

Institute for Infectious Diseases - University of Bern, Switzerland

- Odette J. Bernasconi (PhD student)*
- Fernanda Pimentel (PhD student)
- Thomas Büdel (Lab Technician)
- Dr. med. Baharak Babouee-Flury (PostDoc)



* Hans Sigrist Foundation (Bern, Switzerland)

Kompetenzzentrum für Militär- und Katastrophenmedizin (Switzerland)